



GE Medical Systems

gemedical.com

Technical Publications

Direction 2173225–100

Revision 7

AMX–4+ Service (Model 2169360, 2236420 & 2275938 Series)

Copyright© 1996, 1997, 1999, 2000, 2002 and 2004 By General Electric Co.
All Rights Reserved



THIS PAGE INTENTIONALLY LEFT BLANK.

WARNING

- THIS SERVICE MANUAL IS AVAILABLE IN ENGLISH ONLY.
- IF A CUSTOMER'S SERVICE PROVIDER REQUIRES A LANGUAGE OTHER THAN ENGLISH, IT IS THE CUSTOMER'S RESPONSIBILITY TO PROVIDE TRANSLATION SERVICES.
- DO NOT ATTEMPT TO SERVICE THE EQUIPMENT UNLESS THIS SERVICE MANUAL HAS BEEN CONSULTED AND IS UNDERSTOOD.
- FAILURE TO HEED THIS WARNING MAY RESULT IN INJURY TO THE SERVICE PROVIDER, OPERATOR OR PATIENT FROM ELECTRIC SHOCK, MECHANICAL OR OTHER HAZARDS.

AVERTISSEMENT

- CE MANUEL DE MAINTENANCE N'EST DISPONIBLE QU'EN ANGLAIS.
- SI LE TECHNICIEN DU CLIENT A BESOIN DE CE MANUEL DANS UNE AUTRE LANGUE QUE L'ANGLAIS, C'EST AU CLIENT QU'IL INCOMBE DE LE FAIRE TRADUIRE.
- NE PAS TENTER D'INTERVENTION SUR LES ÉQUIPEMENTS TANT QUE LE MANUEL SERVICE N'A PAS ÉTÉ CONSULTÉ ET COMPRIS.
- LE NON-RESPECT DE CET AVERTISSEMENT PEUT ENTRAÎNER CHEZ LE TECHNICIEN, L'OPÉRATEUR OU LE PATIENT DES BLESSURES DUES À DES DANGERS ÉLECTRIQUES, MÉCANIQUES OU AUTRES.

WARNUNG

- DIESES KUNDENDIENST-HANDBUCH EXISTIERT NUR IN ENGLISCHER SPRACHE.
- FALLS EIN FREMDER KUNDENDIENST EINE ANDERE SPRACHE BENÖTIGT, IST ES AUFGABE DES KUNDEN FÜR EINE ENTSPRECHENDE ÜBERSETZUNG ZU SORGEN.
- VERSUCHEN SIE NICHT, DAS GERÄT ZU REPARIEREN, BEVOR DIESES KUNDENDIENST-HANDBUCH ZU RATE GEZOGEN UND VERSTANDEN WURDE.
- WIRD DIESE WARNUNG NICHT BEACHTET, SO KANN ES ZU VERLETZUNGEN DES KUNDENDIENSTTECHNIKERS, DES BEDIENERS ODER DES PATIENTEN DURCH ELEKTRISCHE SCHLÄGE, MECHANISCHE ODER SONSTIGE GEFAHREN KOMMEN.

AVISO

- ESTE MANUAL DE SERVICIO SÓLO EXISTE EN INGLÉS.
- SI ALGÚN PROVEEDOR DE SERVICIOS AJENO A GEMS SOLICITA UN IDIOMA QUE NO SEA EL INGLÉS, ES RESPONSABILIDAD DEL CLIENTE OFRECER UN SERVICIO DE TRADUCCIÓN.
- NO SE DEBERÁ DAR SERVICIO TÉCNICO AL EQUIPO, SIN HABER CONSULTADO Y COMPRENDIDO ESTE MANUAL DE SERVICIO.
- LA NO OBSERVANCIA DEL PRESENTE AVISO PUEDE DAR LUGAR A QUE EL PROVEEDOR DE SERVICIOS, EL OPERADOR O EL PACIENTE SUFRAN LESIONES PROVOCADAS POR CAUSAS ELÉCTRICAS, MECÁNICAS O DE OTRA NATURALEZA.

ATENÇÃO

- ESTE MANUAL DE ASSISTÊNCIA TÉCNICA SÓ SE ENCONTRA DISPONÍVEL EM INGLÊS.
- SE QUALQUER OUTRO SERVIÇO DE ASSISTÊNCIA TÉCNICA, QUE NÃO A GEMS, SOLICITAR ESTES MANUAIS NOUTRO IDIOMA, É DA RESPONSABILIDADE DO CLIENTE FORNECER OS SERVIÇOS DE TRADUÇÃO.
- NÃO TENHA TENTADO REPARAR O EQUIPAMENTO SEM TER CONSULTADO E COMPREENDIDO ESTE MANUAL DE ASSISTÊNCIA TÉCNICA.
- O NÃO CUMPRIMENTO DESTA AVISO PODE POR EM PERIGO A SEGURANÇA DO TÉCNICO, OPERADOR OU PACIENTE DEVIDO A CHOQUES ELÉTRICOS, MECÂNICOS OU OUTROS.

AVVERTENZA

- IL PRESENTE MANUALE DI MANUTENZIONE È DISPONIBILE SOLTANTO IN INGLESE.
- SE UN ADDETTO ALLA MANUTENZIONE ESTERNO ALLA GEMS RICHIEDE IL MANUALE IN UNA LINGUA DIVERSA, IL CLIENTE È TENUTO A PROVVEDERE DIRETTAMENTE ALLA TRADUZIONE.
- SI PROCEDA ALLA MANUTENZIONE DELL'APPARECCHIATURA SOLO DOPO AVER CONSULTATO IL PRESENTE MANUALE ED AVERNE COMPRESO IL CONTENUTO.
- NON TENERE CONTO DELLA PRESENTE AVVERTENZA POTREBBE FAR COMPIERE OPERAZIONI DA CUI DERIVINO LESIONI ALL'ADDETTO ALLA MANUTENZIONE, ALL'UTILIZZATORE ED AL PAZIENTE PER FOLGORAZIONE ELETTRICA, PER URTI MECCANICI OD ALTRI RISCHI.

警告

- ・このサービスマニュアルには英語版しかありません。
- ・GEMS以外でサービスを担当される業者が英語以外の言語を要求される場合、翻訳作業はその業者の責任で行うものとさせていただきます。
- ・このサービスマニュアルを熟読し理解せずに、装置のサービスを行わないで下さい。
- ・この警告に従わない場合、サービスを担当される方、操作員あるいは患者さんが、感電や機械的又はその他の危険により負傷する可能性があります。

注意:

- 本维修手册仅存有英文本。
- 非 GEMS 公司的维修员要求非英文本的维修手册时，客户需自行负责翻译。
- 未详细阅读和完全了解本手册之前，不得进行维修。
- 忽略本注意事项会对维修员，操作员或病人造成触电，机械伤害或其他伤害。

IMPORTANT! . . . X-RAY PROTECTION



X-ray equipment if not properly used may cause injury. Accordingly, the instructions herein contained should be thoroughly read and understood by everyone who will use the equipment before you attempt to place this equipment in operation. The General Electric Company, Medical Systems Group, will be glad to assist and cooperate in placing this equipment in use.

Although this apparatus incorporates a high degree of protection against x-radiation other than the useful beam, no practical design of equipment can provide complete protection. Nor can any practical design compel the operator to take adequate precautions to prevent the possibility of any persons carelessly exposing themselves or others to radiation.

It is important that everyone having anything to do with x-radiation be properly trained and fully acquainted with the recommendations of the National Council on Radiation Protection and Measurements as published in NCRP Reports available from NCRP Publications, 7910 Woodmont Avenue, Room 1016, Bethesda, Maryland 20814, and of the International Commission on Radiation Protection, and take adequate steps to protect against injury.

The equipment is sold with the understanding that the General Electric Company, Medical Systems Group, its agents, and representatives have no responsibility for injury or damage which may result from improper use of the equipment.

Various protective material and devices are available. It is urged that such materials or devices be used.

CAUTION: United States Federal law restricts this device to use by or on the order of a physician.

CERTIFIED ELECTRICAL CONTRACTOR STATEMENT



All electrical installations that are preliminary to positioning of the equipment at the site prepared for the equipment shall be performed by licensed electrical contractors. In addition, electrical feeds into the Power Distribution Unit shall be performed by licensed electrical contractors. Other connections between pieces of electrical equipment, calibrations, and testing shall be performed by qualified GE Medical personnel. The products involved (and the accompanying electrical installations) are highly sophisticated, and special engineering competence is required. In performing all electrical work on these products, GE will use its own specially trained field engineers. All of GE's electrical work on these products will comply with the requirements of the applicable electrical codes.

The purchaser of GE equipment shall only utilize qualified personnel (i.e., GE's field engineers, personnel of third-party service companies with equivalent training, or licensed electricians) to perform electrical servicing on the equipment.

DAMAGE IN TRANSPORTATION

All packages should be closely examined at time of delivery. If damage is apparent, have notation "**damage in shipment**" written on **all** copies of the freight or express bill **before** delivery is accepted or "signed for" by a General Electric representative or a hospital receiving agent. Whether noted or concealed, damage **MUST** be reported to the carrier **immediately** upon discovery, or in any event, within **14** days after receipt, and the contents and containers held for inspection by the carrier. A transportation company will not pay a claim for damage if an inspection is not requested within this **14** day period.

Call Traffic and Transportation, Milwaukee, WI (414) 827-3468 / 8*285-3468 **immediately** after damage is found. At this time be ready to supply name of carrier, delivery date, consignee name, freight or express bill number, item damaged and extent of damage.

Complete instructions regarding claim procedure are found in Section "S" of the Policy & Procedure Bulletins.

6/17/94

THIS PAGE INTENTIONALLY LEFT BLANK.

REVISION HISTORY

REV	DATE	REASON FOR CHANGE
A	Nov. 6 1996	Initial Draft.
0	Dec. 13, 1996	Initial release.
1	Mar. 7, 1997	Made cable drape illustration changes, brake relay illustration changes; updated battery installation illustration.
2	June 24, 1997	Modified vertical arm latch switch to new style, added latch wiring diagram in Section 9.
3	Aug. 14, 1997	High Impact Inspection.
4	Apr. 12, 1999	Added AMX-4+ Model 2236420 Series.
5	Nov. 8, 2000	Added AMX-4+ Model 2275938 Series.
6	Sept. 10, 2002	Added NOTE to section 4-2 about Loctite and cleaning screw holes.
7	Oct. 21, 2004	Added model 46-270615P3 collimator information to replacement procedures.

THIS PAGE INTENTIONALLY LEFT BLANK.

TABLE OF CONTENTS

SECTION 1 – SERVICING THE COLLIMATOR	1
1-1 Collimator Removal	1
1-2 Collimator Interface and Port Plates Check	2
1-3 Collimator Installation	4
1-4 Tape Measure Replacement	4
1-5 Field Lamp Replacement	5
1-6 Light Field Adjustment	5
1-6-1 Adjustment Method 1	6
1-6-2 Adjustment Method 2	8
1-7 Crosshair Window Adjustment	10
1-8 Field Light Switch Replacement	10
1-9 Field Size Pointer Replacement	10
1-10 Front Panel Replacement (Medys collimator)	11
1-11 Brake Switch Replacement	12
SECTION 2 – SERVICING THE X-RAY TUBE	14
2-1 X-Ray Tube Removal	14
2-2 X-Ray Tube Installation	15
SECTION 3 – SERVICING THE HORIZONTAL ARM	16
3-1 Yoke Assembly Removal	16
3-2 Yoke Assembly Installation	17
3-3 Yoke Mounting Plate Assembly Removal	17
3-4 Yoke Mounting Plate Assembly Installation	19
3-5 Horizontal Arm Assembly Removal	20
3-6 Horizontal Arm Assembly Installation	21
3-7 Horizontal Arm Brake and Cable Assembly Removal	22
3-8 Horizontal Arm Brake and Cable Assembly Installation	22
3-9 Horizontal Arm Brake Cable Removal	24
3-10 Horizontal Arm Brake Cable Installation	24
3-11 Horizontal Arm Brake Removal	25
3-12 Horizontal Arm Brake Installation	25
SECTION 4 – SERVICING THE VERTICAL COLUMN	26
4-1 Remove Vertical Column Assembly	26
4-2 Install Vertical Column Assembly	27
4-3 Counterweight Cables	28
4-3-1 Main Cable	29
4-3-2 Safety Cable	31
4-3-3 Reassemble Unit	31
4-3-4 Vertical Column Cable Tension	31
4-3-5 Remove Follower Cable	32
4-3-6 Install Follower Cable	34
4-4 Upper Brake, Pulley, and Bearings	34
4-4-1 Disconnect Brake	34
4-4-2 Remove Upper Pulley and Brake	34
4-4-3 Replace Upper Pulley and Bearings	35
4-4-4 Install Upper Pulley, Bearings and Brake	36
4-4-5 Connect Brake	37
4-4-6 Remove Upper Pulley, Bearings and Brake (Alternate)	37
4-4-7 Install Upper Pulley, Bearings and Brake (Alternate)	37

TABLE OF CONTENTS – Continued

4-5	Counterweight	38
4-5-1	Remove Counterweight	38
4-5-2	Install Counterweight	39
4-5-3	Counterweight Balancing	39
4-6	Carriage Assembly	40
4-6-1	Remove Carriage	40
4-6-2	Install Carriage	40
4-6-3	Remove Carriage Bearings	40
4-6-4	Install Carriage Bearings	41
4-7	Lower Brake Removal	42
4-8	Lower Bearings and Baseplate Removal	42
4-9	Lower Brake and Bearings Installation	43
4-10	Column Base Removal	43
4-11	Column Base Installation	43
4-12	Lower Pulley and Bearings Removal	44
4-13	Lower Pulley and Bearings Installation	44
4-14	Column Support Assembly Removal	44
4-15	Column Support Assembly Installation	45

SECTION 5 – SERVICING THE CASSETTE DRAWER, COVERS, HORIZONTAL ARM LATCH AND HANDSWITCH

5-1	Cassette Drawer Removal	46
5-2	Cassette Drawer Installation	46
5-3	Side Covers Removal	47
5-4	Side Cover Installation	47
5-5	Top Cover Removal	48
5-6	Top Cover Installation	48
5-7	Remove Latch Assembly and Components (Hall Effect Switch)	49
5-8	Install Latch Assembly and Components (Hall Effect Switch)	51
5-9	Remove And Install Hand Switch Cable	52
5-10	Install Hand Switch	52

SECTION 6 – SERVICING THE BATTERIES

6-1	Remove Battery Cover	54
6-2	Disconnect Battery Sense Harness (If Required)	55
6-3	Disconnect and Remove Batteries	55
6-4	Disconnect and Remove Battery Test Terminal Strip	58
6-5	Continue Battery Disconnection/Removal	58
6-6	Install Batteries	58
6-7	Prepare to Install Batteries	59
6-8	Install and Connect Batteries	59
6-9	Install and Connect Battery Test Terminal Strip	61
6-10	Continue Battery Installation/Connection	61
6-11	Connect Battery Terminal Test Strip	61
6-12	Continue Battery Installation/Connection	61
6-13	Connect Battery Terminal Test Strip	61
6-14	Complete Battery Installation/Connection	62
6-15	Install Battery Compartment Cover	62

TABLE OF CONTENTS – Continued

SECTION 7 – SERVICING THE BASE ASSEMBLY	63
7-1 Front Bumper	63
7-1-1 Remove Pads	63
7-1-2 Install Pads	64
7-1-3 Bumper Switch Removal	64
7-1-4 Bumper Switch Installation	64
7-1-5 Bumper Bearings Removal	65
7-1-6 Bumper Bearings Installation	65
7-2 Raise Unit for Service	66
7-2-1 Raise Rear Wheels	66
7-2-2 Lower Rear Wheels	66
7-2-3 Raise Front Wheels	67
7-2-4 Lower Front Wheels	68
7-3 Remove Drive Wheels	68
7-4 Install Drive Wheels	69
7-5 Remove Motor and Brake	70
7-6 Install Motor and Brake	71
7-7 Remove Caster Wheel	72
7-8 Install Caster Wheel	72
SECTION 8 – SERVICING THE DRIVE HANDLE ASSEMBLY	73
8-1 Handle Sensor Removal	73
8-2 Handle Sensor Installation	74
8-3 Handle Drive Enable Switch Removal	75
8-4 Handle Drive Enable Switch Installation	75
SECTION 9 – SERVICING THE HIGH VOLTAGE TRANSFORMER	77
9-1 High Voltage Transformer Removal	77
9-2 High Voltage Transformer Installation	78
SECTION 10 – SERVICING THE POWER CORD AND REEL ASSEMBLY AND CABLE ROUTING .	79
10-1 Cord/Cord Reel Removal	79
10-1-1 Cord Reel Removal (Domestic Models)	79
10-1-2 Cord Reel Removal (IEC Models)	79
10-2 Cord/Cord Reel Installation	80
10-2-1 Cord Reel Installation (Domestic Models)	80
10-2-2 Cord Reel Installation (IEC Models)	80
10-3 Horizontal Arm and Vertical Column Cable Routing	80
10-3-1 Cathode Cable	81
10-3-2 Anode Cable	83
10-3-3 Collimator Cable	84
10-3-4 Stator Cable	88
10-3-5 Horizontal Arm Brake Cable	89
10-3-6 Control Cable	89
10-3-7 Vertical Column Brake Cable	90

THIS PAGE INTENTIONALLY LEFT BLANK.


AMX-4+ SERVICE**Identification**

See Illustration 9. The AMX-4+ is identified by Model Number on the rating plate located on the top cover. Model part and catalog numbers are identified in Table 1. This direction contains instructions for performing mechanical maintenance on these Units.

TABLE 1
AMX-4+ MODELS

Description	Part number	catalog number	Part number	catalog number
domestic	2169360-7	a0659f	2236420-7 & 2275938-7	A0659JF
domestic, aec	2169360-8	a0659fa	2236420-8 & 2275938-8	A0659JG
domestic, tech switch	2169360-9	a0659fc	2236420-9 & 2275938-9	A0659JH
domestic, aec, tech switch	2169360-10	a0659fb	2236420-10 & 2275938-10	A0659JJ
iec, emc	2169360	a0659a	2236420 & 2275938	A0659J
iec, emc, aec	2169360-2	a0659aa	2236420-2 & 2275938-2	A0659JA
iec, emc, tech switch	2169360-3	a0659ab	2236420-3 & 2275938-3	A0659JB
iec, emc, aec, tech switch	2169360-4	a0659ac	2236420-4 & 2275938-4	A0659JC
japan	2169360-5	a0659c	2236420-5 & 2275938-5	A0659JD
japan short column	2169360-6	a0659d	2236420-6 & 2275938-6	A0659JE

ILLUSTRATION 9
AMX-4+ IDENTIFICATION



The instructions cover the removal, installation, and adjustment of the field serviceable assemblies and parts comprising the unit. Illustrations are provided as a supplement to the procedures to assist the service technician in understanding the relationship of the parts involved.

Definition of Terms

The front end of the unit is defined as the end containing the vertical column and bumper assembly. Left and right sides are determined by standing at the rear (handle end) of the unit, facing the unit.

Tools and Materials

A list of tools required to perform the service adjustment procedures is given below. With the exception of the vertical column lifting tool and floor jack, the tools listed are part of the standard complement of tools supplied in the service technician's tool kit.

Tools

- Vertical column lifting tool (46-303362P1)
- AMX-4 Series Lifting Jack (46-302966P1)
- Counterweight Lifting Tool (46-303491G1)
- HHS Field Test Kit (46-177372G1)
- Torque Wrench 2-36 pound-inches/0.23 4.1 Newton-meters (46-282782P1)
- Torque Wrench Kit (46-268445G1)
- H/V Cable Spanner Wrench (507A935G1)

Materials

- Loctite sealant type 242 (46-170686P2 10 cc bottle) and 290 (46-170687P2 50cc bottle)
- Bearing grease - Lithium (Fiske 630-A-A) or equal
- Silicone insulating compound (46-125224P3 3 oz. tube)
- Diala AX Dielectric Insulating Oil (T0552G 1 gallon)
- Electrical tape
- Cotton rags
- Cotton swabs

Periodic Inspection

Periodic inspection instructions for various assemblies of the AMX-4+ are listed in the maintenance schedule given in Direction 2173227-100, *AMX-4+ Periodic Maintenance*. All deficiencies uncovered during inspection should be corrected before further use of the equipment.

Removal and Installation Procedure Index

A list of AMX-4+ assemblies for which removal and installation procedures are provided in the following index. The index lists the section in which the procedure is presented, the name of the assembly, and the number of the paragraph that contains the removal or installation procedure.

SECTION	ASSEMBLY NAME	REMOVAL PARAGRAPH	INSTALLATION PARAGRAPH
TITLE			
Section 1 Servicing the Collimator Assembly	Collimator Assembly	1-1	1-2
Section 2 Servicing the X-ray Tube Assembly	X-Ray Tube Assembly	2-1	2-2
Section 3 Servicing the Horizontal Arm Assembly	Yoke Assembly Yoke Mounting Plate Assembly Horizontal Arm Assembly Brake – Cable Assembly Cable Assembly Brake Assembly	3-1 3-3 3-5 3-7 3-9 3-11	3-2 3-4 3-6 3-8 3-10 3-12
Section 4 Servicing the Vertical Column Assembly	Vertical Column Cables Upper Pulley and Brake Assembly Counterweights Upper Brake Assembly Upper Pulley and Bearings Carriage Assembly Carriage Bearings Lower Brake and Bearings Column Base Lower Pulley and Bearings Column Support	4-1 4-3 4-5 4-7 4-9 4-11 4-13 4-15 4-17 4-19 4-21 4-23	4-2 4-4 4-6 4-8 4-10 4-12 4-14 4-16 4-18 4-20 4-22 4-24
Section 5 Servicing the Trim Covers	Cassette Drawer Side Covers Top Cover	5-1 5-3 5-5	5-2 5-4 5-6
Section 6 Servicing the Base Assembly	Gates Battery Bumper Assembly Bumper Switch Bumper Bearings Motor and Brake Drive Wheel Caster Wheel	6-1 6-3 6-5 6-7 6-9 6-11 6-13	6-2 6-4 6-6 6-8 6-10 6-12 6-14
Section 7 Servicing the Handle Assembly	Sensor Assembly Drive Enable Switch	7-1 7-3	7-2 7-4
Section 8 Servicing the High Voltage Transformer	Transformer	8-1	8-2
Section 9 Servicing the Latch Assembly and Cable Routing	Latch Assembly (Hall Effect Switch) Cord/Cord Reel Cord Reel Cable Routing Cathode Cable Anode Cable Collimator Cable Stator Cable Horizontal Arm Brake Cable Control Cable Vertical Column Brake Cable	9-1 9-3 9-3-2	9-2 9-4 9-5 9-5-1 9-5-2 9-5-3 9-5-4 9-5-5 9-5-6 9-5-7

THIS PAGE INTENTIONALLY LEFT BLANK.

SECTION 1 – SERVICING THE COLLIMATOR

This section provides procedures covering the removal and installation of the collimator. It also provides procedures for replacing and aligning the field lamp, aligning the crosshairs, and replacing the brake and field lamp switches.

1-1 Collimator Removal


1. Close collimator blades.
2. Rotate X-ray tube and collimator assembly so collimator crosshair window is facing upward.
3. Place the horizontal arm assembly in park position on top cover of unit.

CAUTION


Potential for sudden and unexpected arm motion: Secure the horizontal arm before removing components from the arm. Removing weight allows the counterweight to descend rapidly and the arm to ascend rapidly possibly causing damage to the unit and personal injury. Use Counterweight Lift Tool 46-303491G1.

4. Shut off power by turning the main circuit breaker to the OFF position.
5. Remove two screws in bottom of terminal box located on back of vertical column and lift off cover. See Illustration 1-2.
6. Disconnect collimator cable leads from terminal block. Record and tag leads.
7. Remove Heyco bushing and cable through bottom of terminal box.

ILLUSTRATION 1-2
TERMINAL BOX AND COLLIMATOR CABLE ROUTING



8. Cut ty-raps securing collimator cable to cathode cable and remove from Strap Points G and N, and Clamp Points E and H.
9. Remove three 10-24 $\times$ 0.75 hex socket head cap screws securing collimator assembly to X-ray tube. See Illustration 1-3.
10. Remove collimator from interface plate.

ILLUSTRATION 1-3
COLLIMATOR

1-2 Collimator Interface and Port Plates Check

Note: The collimator interface plate for every new and reloaded tube from the GE Reload Facility is factory aligned and should not be disturbed. If plate alignment is altered, the plate will have to be re-aligned in the field.

1. Close collimator blades.
2. Rotate X-ray tube and collimator assembly so collimator crosshair window is facing upward.
3. Place the horizontal arm assembly in park position on top cover of unit.

CAUTION


Potential for sudden and unexpected arm motion: Secure the horizontal arm before removing components from the arm. Removing weight allows the counterweight to descend rapidly and the arm to ascend rapidly possibly causing damage to the unit and personal injury. Use Counterweight Lift Tool 46-303491G1.

4. Shut off power by turning the main circuit breaker to the OFF position.
5. Remove three 10-24 $\times$ 0.75 inch hex socket head cap screws securing collimator assembly to X-ray tube. See Illustration 1-3.

**ILLUSTRATION 1-4
PORT AND COLLIMATOR INTERFACE PLATE INSTALLATION**

MOBILE SYSTEMS USE EITHER
FOUR OR EIGHT 1/4-20 X 0.625 INCH
PORT PLATE SCREWS

10-32 X 0.50 INCH
INTERFACE PLATE
SCREWS



6. Remove collimator from interface plate.
7. Carefully place collimator and cables on horizontal arm.
8. Check that each X-ray port plate screw is securely fastened in place.
9. If any port plate screw is loose, unscrew it until it touches the interface plate.

NOTICE

Do not remove collimator interface plate. This plate is factory aligned and should not be disturbed. (However, if any of the interface plate screws are loose, re-torque them to 32–40 pound–inches (3.6–4.5 N–m)).

10. Apply Loctite 290 to the port plate screw and tighten to 45–50 pound–inches (5.1– 5.6 N–m). Loctite 290 is thin and, like a wick, will draw down around the screw threads. Clean the area of any residue solution.
11. Turn collimator upside down and insert cone into the tube port interface plate. Align mounting screw holes.
12. Attach collimator assembly to tube port interface plate with three 10–24 × 0.75 inch hex socket head cap screws using Loctite 242. Torque to 50–60 pound–inches (5.6–6.8 N–m).

1-3 Collimator Installation

1. Disconnect power from the unit by switching the main circuit breaker to the OFF position.
2. Horizontal arm assembly is in parked position on top cover of unit (or counterweight is secured)
3. X–ray tube is facing upward.

NOTICE

Equipment Damage: Do not bump the cone when handling the collimator, especially when mounting it on the tube. It is made of lead and damages easily.

4. Turn collimator upside down and insert cone into the tube port interface plate. Align mounting screw holes.
5. Attach collimator assembly to tube port interface plate with three 10–24 × 0.75 inch hex socket head cap screws using Loctite 242. Torque to 50–60 pound–inches (5.6–6.8 N–m).
6. Route collimator cable along the cathode cable to strap point on side of the column. Secure collimator cable in place with ty–raps to the cathode cable. See Illustration 1–2.
7. Drape the cable from the side of the column to the terminal box on the back of the column. Insert the cable with a Heyco bushing up through the bottom of the terminal box.
8. Connect collimator cable leads to terminal block in accordance with connection data recorded and tagged on leads.
9. Position back cover on terminal box and secure with two screws.
10. Rotate X–ray tube and collimator assembly to a port down position.

Note: If the collimator being installed is a replacement, perform “Beam Quality (Half Value Layer)” in Direction 2173227–100 *AMX–4+ Periodic Maintenance*. Also perform “Functional Test – Manual Rad Collimator Version” in Tab 6 of Direction 46–013894, *System Field Test for HHS*.

1-4 Tape Measure Replacement

Note: This procedure is for the **Micro Medical Systems** collimator (46–270615P2) only. For the **Medys** collimator (46–270615P3), the tape measure is not a replaceable part.

1. Remove the two (2) Phillips screws (M3x6) which attach the old tape measure to the left–side collimator cover.
2. Install the new tape measure and tighten the screws:
 - For a HRT X–ray tube (46–155750G804), use the upper mounting holes of the tape measure.
 - For an Orion X–ray tube (2185226), use the lower mounting holes of the tape measure.

1-5 Field Lamp Replacement


1. Disconnect power from the unit by switching the main circuit breaker to the OFF position.
2. Remove the rear cover on the collimator by removing four metric Phillips head screws holding it in place. Remove the cover. See Illustration 1-5.
3. Close the collimator blades. This makes it easier to reach the lamp for removal and replacement.

NOTICE

Potential for equipment damage (field lamp): Field lamp should never be touched with bare fingers since chemicals on skin will damage lamp and reduce its useful life. Always use gloves, tissue or soft cloth when handling lamp. If the outside of the lamp becomes contaminated, clean with a soft cloth dampened in water. Wipe dry after cleaning.

4. Carefully remove the lamp by pulling it straight back and out of its receptacle. Replace it with a new lamp. Refer to renewal parts listing for correct lamp.
5. Whenever a lamp is replaced, the light field intensity must be checked. Refer to the HHS Collimator Test section of Direction 46-013894 (System Field Test for HHS).
6. Whenever a lamp is replaced, the light field and X-ray field sizes must be checked to be certain the difference between them does not exceed a specified maximum tolerance. Refer to Section 1-6, *Light Field Adjustment*, to make these measurements and any necessary adjustments.

ILLUSTRATION 1-5
FIELD LAMP REPLACEMENT



1-6 Light Field Adjustment

The X-ray field and light field must be properly aligned to meet HHS requirements. Two methods for accomplishing this are given below.


Method 1, Section 1-6-1 is only used if the "Light Field to X-Ray Field Test Pattern," 46-198466P1, from an HHS Field Test Kit, 46-177372G1, is available.

Method 2, Section 1-6-2, uses lead strips to indicate the light field areas.

1-6-1 Adjustment Method 1

1. Place loaded cassette (14×17 inches suggested) on floor directly beneath collimator.
2. Adjust the collimator so that it is parallel to the cassette.
3. Set SID to 40 inches.
4. Place light field to X-ray field test pattern on top of cassette.
5. Turn on field lamp.
6. Adjust collimator blades so light field exactly covers the 8×10 inch (20.32×25.4 cm) rectangle on the test pattern. See Illustration 1-6.

ILLUSTRATION 1-6

APPEARANCE OF X-RAY FIELD AND LIGHT FIELD ALIGNMENT ON EXPOSED FILM USING TEST PATTERN

7. Expose film.
 - > 60 kVp, 2 mAs.
 - > Do not disturb the cassette and pattern because a second exposure is to be made.
8. Open collimator blades wide enough to include test pattern.
9. Expose film a second time and develop the film.
 - > 60 kVp, 2 mAs.
 - > This exposure should include the 8×10 inch rectangle. Refer to Illustration 1-4. If the rectangle does not show, reload the cassette and repeat steps 1 through 8 but open the collimator blades in step 7 wider than last time.
10. Measure and record the horizontal A & B distances, in inches, between the X-ray field and the light field. See Illustration 1-6.
11. Measure and record the vertical C & D distances, in inches, between the X-ray field and the light field. See Illustration 1-6.

12. The total of the horizontal differences as well as the total of the vertical differences must each be less than 1.8% of the SID:

$$> A + B < 1.8\% \text{ SID} = A + B < 0.018 \times 40 \text{ inches} = \underline{A + B < 0.72 \text{ inches}}$$


$$> C + D < 1.8\% \text{ SID} = C + D < 0.018 \times 40 \text{ inches} = \underline{C + D < 0.72 \text{ inches}}$$

If the above requirements are met, go to step 18. If they are not met, proceed with the following steps.

Note: Ideally, the X-ray field and the light field should coincide but this is rarely the case. The best compromise is to move the light field so the X-ray field is centered within it. Usually the light field area will be larger than the X-ray area. The size and location of the light field can be adjusted by moving the bulb.

13. To place the two fields as close to coincidence as possible, only the light field area is moved. The light field area adjustment is made by repositioning the field lamp housing.
14. Loosen the two Phillips head adjustment screws located in the back of the collimator. See Illustration 1-7.

ILLUSTRATION 1-7
FIELD LAMP ADJUSTMENT



15. Examine the last exposure and determine which direction to move the field lamp to best position it around the X-ray field.
16. Turn on the field lamp.

CAUTION

Bright light might cause temporary blindness: Components are hot in the high energy field lamp adjustment area. Do not look directly at the field lamp.

17. The field lamp adjustment area will be hot, so use a metal tool to move the bracket when repositioning the light field area on the film cassette.
18. When the light field area has been repositioned equal distant around the X-ray field, tighten the adjustment screws. Alternately tightening the two screws will help prevent the bracket from slipping out of place.
19. Repeat steps 1 through 16 until the "percentage of SID" requirements are obtained.

Note: If the percentage of SID requirement cannot be obtained, the light field is either too large or in rare instances too small. In this case the collimator will have to be returned to the factory.


20. Replace the back cover and fasten it in place with four metric Phillips head screws. See Illustration 1-5.
21. After any adjustment to the field light area, the collimator crosshair window position should be checked and corrected if necessary. Refer to Section 1-5.

1-6-2 Adjustment Method 2

1. Place loaded cassette (14 × 17 inches (35.56 x 43.18 cm) suggested) on floor directly beneath collimator.
2. Set SID to 40 inches.
3. Set collimator for 8 × 10 inches.
4. Turn on field lamp.
5. Place lead strips at each corner of light field to indicate size of field. See Illustration 1-8.
6. Expose film. Do not move film cassette.
 - > 60 kVp, 2 mAs.
7. Open collimator blades.
8. Expose film a second time and develop the film.
 - > 60 kVp, 2 mAs.

Note: This exposure should include the lead strips marking the corners of the light field. See Illustration 1-8. If the lead strips do not show, reload the cassette and repeat steps 1 through 8 but open the collimator blades in step 7 far enough to include the lead strips.

ILLUSTRATION 1-8
APPEARANCE OF X-RAY FIELD AND LIGHT FIELD ALIGNMENT ON EXPOSED FILM USING LEAD STRIPS



9. Measure and record the horizontal A & B distances, in inches, between the X-ray field and the light field. See Illustration 1-8.
10. Measure and record the vertical C & D distances, in inches, between the X-ray field and the light field. See Illustration 1-7.

11. The total of the horizontal differences as well as the total of the vertical differences must each be less than 1.8% of the SID:

$$> A + B < 1.8\% \text{ SID} = A + B < 0.018 \times 40 \text{ inches} = \underline{A + B < 0.72 \text{ inches}}$$


$$> C + D < 1.8\% \text{ SID} = C + D < 0.018 \times 40 \text{ inches} = \underline{C + D < 0.72 \text{ inches}}$$

If the above requirements are met, go to step 18. If they are not met, proceed with the following steps.

Note: Ideally, the X-ray field and the light field should coincide but this is rarely the case. The best compromise is to move the light field so the X-ray field is centered within it. Usually the light field area will be larger than the X-ray area.

12. To place the two fields as close to coincidence as possible, only the light field area is moved. The light field area adjustment is made by repositioning the field lamp housing.
13. Loosen the two Phillips head adjustment screws located in the back of the collimator. See Illustration 1-9.
14. Examine the last exposure and determine which direction to move the field lamp to best position light field around the X-ray field.
15. Turn on the field lamp.

ILLUSTRATION 1-9
FIELD LAMP ADJUSTMENT



CAUTION

Bright light might cause temporary blindness: Components are hot in the high energy field lamp adjustment area. Do not look directly at the field lamp.

16. The field lamp adjustment area will be hot so use a metal tool to move the bracket when repositioning the light field area on the film cassette.
17. When the light field area has been repositioned equidistant around the X-ray field, tighten the adjustment screws. Alternately tightening the two screws will help prevent the bracket from slipping out of place.
18. Repeat steps 1 through 16 until the percentage of SID requirements are obtained.

Note: If the percentage of SID requirement cannot be obtained, the light field is either too large or in rare instances too small. In this case the collimator will have to be returned to the factory.

19. Replace the back cover and fasten it in place with four metric Phillips head screws. See Illustration 1-8.
20. After any adjustment to the field light area, the collimator crosshair window position should be checked and corrected if necessary. Refer to Section 1-7.

1-7 Crosshair Window Adjustment

1. After any adjustment to the field light area the collimator crosshair window position should be checked and corrected if necessary.
2. Place Light Field to X-Ray Field Test Pattern, 46-198466P1, directly under the collimator. If the test pattern is not available use a sheet of rectangular paper and draw a cross extending from exact center to each edge and parallel to the sides of the sheet.
3. Turn on field lamp and adjust collimator blades so edges of light field are coincident with the rectangle on the test pattern or the sides of the paper sheet. This may involve moving the test pattern or paper sheet.
4. Check the window crosshairs. The crosshairs must be within 1/4 inch of the pattern centerlines or cross (including the width of the crosshair shadow) over the entire length of the crosshairs.
5. If the window crosshairs are not properly aligned, go to step 7. If the window crosshairs are properly aligned, do not perform steps 7 through 9.
6. Loosen the four screws retaining the crosshair window in the bottom or exit side of the collimator.
7. Gently slide the window until the crosshairs are aligned with the test pattern centerlines or the cross on the paper.
8. Tighten the four screws retaining the crosshair window.

1-8 Field Light Switch Replacement

Note: This procedure is for the **Micro Medical Systems** collimator (46-270615P2) only. For the **Medys** collimator (46-270615P3), the field light switch is not a replaceable part.

1. Place horizontal arm in lock position.
2. Shut off power by turning the main circuit breaker to the OFF position.
3. Use a small screwdriver and very carefully pry the field light switch out of the collimator control panel. The switch press fits into the panel and will easily snap out of the panel. A short length of two wires will be attached to the switch. Use care not to scratch the panel during this procedure. See Illustration 1-3.
4. Unsolder the two wires.
5. Resolder two wires to the terminals on new replacement switch.
6. Carefully insert the switch and wires into the panel and press into place.

1-9 Field Size Pointer Replacement

Note: This procedure is for the **Micro Medical Systems** collimator (46-270615P2) only. For the **Medys** collimator (46-270615P3), the field size pointers are replaceable parts.

1. Use a 2 mm hex wrench to loosen the two (2) M4 setscrews that attach the knob to the shaft. Slide the knob off of the shaft.
2. Use a 1.5 mm hex wrench to loosen the two (2) M3 setscrews that attach the pointer to the shaft. Slide the pointer off of the shaft.

3. Slide the new pointer onto the shaft and lightly tighten the two set screws. Perform the following steps to correctly align the pointer:
 - a. Place a loaded 12 inch x 14 inch cassette (or phosphor tool – 2120565) on the floor and center it to the collimator light field.
 - b. Use the tape measure to position the collimator/X-ray tube at a 40 inch SID.
 - c. Adjust the collimator so that the pointers are set for a 10 inch x 10 inch field.



POTENTIAL EXISTS FOR EXPOSURE TO IONIZING RADIATION. USE A LEAD APRON OR STAND BEHIND LEAD SHIELDING DURING EXPOSURE.

- d. Select 60kVp, 4–5 mAs (or follow exposure instructions on phosphor tool) and make an exposure.
- e. Measure the exposed area size on the developed film (or read from the phosphor tool).
- f. Verify that the indicated field size (as shown by the pointer) matches the actual exposed area size within 1.7% SID (0.68 inches).
- g. If necessary, loosen pointer setscrews, adjust pointer to match the measured size, then tighten the setscrews. Take another exposure to verify that the indicated field size (as shown by the pointer) matches the actual exposed area size within 1.7% SID (0.68 inches). Repeat this step until the specification is met.
- h. Install the knob onto the shaft and tighten the two setscrews.

1-10 Front Panel Replacement (Medys collimator)

Note: This procedure is for the **Micro Medical Systems** collimator (46–270615P2) only. For the **Medys** collimator (46–270615P3), the front panel is not a replaceable part.

1. Shut off power by turning the main circuit breaker to the OFF position.
2. Close the collimator blades.
3. Rotate the X-ray tube and collimator assembly so the collimator cross-hair window is facing upward.
4. Place horizontal arm in lock position.
5. Place the horizontal arm assembly in park position on the top cover of the unit.



Potential for sudden and unexpected arm motion: Secure the horizontal arm before removing components from the arm. Removing weight allows the counterweight to descend rapidly and the arm to ascend rapidly possibly causing damage to the unit and personal injury. Use Counterweight Lift Tool 46–303491G1

6. Disassembly:
 - a. Use a 2 mm hex wrench to loosen the two (2) M4 setscrews that attach the knob to the shaft. Slide the knob off of the shaft.
 - b. Use a 1.5 mm hex wrench to loosen the two (2) M3 setscrews that attach the pointer to the shaft. Slide the pointer off of the shaft.
 - c. Remove the rear cover on the collimator by removing four metric Phillips head screws holding it in place. Remove the rear cover.
 - d. Remove the two (2) Phillips screws (M3x6) which attach the old tape measure to the left-side collimator cover. Remove the tape measure.
 - e. Remove the four (4) screws M3x6 (Phillips) fastening the lateral cover to the top of the central body.

- f. Remove the two (2) screws M3x6 (Phillips) fastening the lateral cover to the rear of the central body.
- g. Remove the lateral cover by pulling it towards the rear of the collimator.
- h. Unscrew the six (6) screws M3x6 (Phillips) fastening the front panel to the left and right sides of the central body, but **do not remove the front panel yet.**


Note: The screws are secured with thread fastener and may require additional effort to unscrew them.

- i. Carefully examine the routing for the wires that are connected to the front panel so that you will be able to duplicate the routing when the new front panel is attached to the central body. If the cable routing is not correct, the front panel is difficult to install.
 - j. Carefully pull out the front panel. **Remember; the front panel is still linked to the central body by means of the switches connector.**
 - k. Disconnect the switches connector and remove the front panel.
7. Install the new front panel (reverse the disassembly sequence). If the new front panel does not line up with the six (6) mounting holes in the central body, make sure that the cable routing is the same as noted in step i (above).
 8. Final Checks – When the new front panel is assembled with its knobs and pointers, perform the following checks and adjustments:
 - a. Light Field adjustment (see Section 1-6 in this manual).
 - b. Cross hair Window adjustment (see Section 1-7 in this manual).
 - c. Field Size Indicator adjustment (see Section 1-9 in this manual).

1-11 Brake Switch Replacement

Note: This procedure is for the **Micro Medical Systems** collimator (46-270615P2) only. For the **Medys** collimator (46-270615P3), the brake switch is not a replaceable part.

1. Place horizontal arm in lock position.
2. Shut off power by turning the main circuit breaker to the OFF position.
3. To prevent parts from falling out of the collimator handles when opening them, rotate the X-ray tube to position the back side of the collimator handles pointing upward.
4. Select the switch to be replaced and remove the two Phillips head screws from the switch retaining plate. Remove the plate.
5. Lift off the switch button and button return spring.
6. Pull the switch from the handle.
7. Unsolder the two wires.
8. Solder the two wires onto the new replacement switch.
9. Place the new switch in the same position as the switch which was replaced.
10. Replace the button return spring and the switch button.
11. Replace the switch retaining plate. Install and tighten the two Phillips head screws.
12. Test the switch to be sure that it activates/deactivates, and does not stick in the ON position.

ILLUSTRATION 1-10
BRAKE SWITCH REPLACEMENT

SECTION 2 – SERVICING THE X-RAY TUBE

This section provides procedures covering the removal and installation of the X-ray tube.


2-1 X-Ray Tube Removal

CAUTION

Secure the horizontal arm before removing components from the arm. Removing weight allows the counterweight to descend rapidly and the arm to ascend rapidly possibly causing damage to the unit and personal injury. Use the Counterweight Lift Tool 46-303491G1.

1. Remove collimator from X-ray tube. Refer to Section 1-1, steps 1 through 10 for removal instructions.
2. Disconnect anode and cathode High Voltage cable connectors from X-ray tube. See Illustration 2-1.

ILLUSTRATION 2-1
X-RAY TUBE ASSEMBLY



3. Remove four button head screws securing top and rear of back cover to X-ray tube and lift off cover.
4. Disconnect control cable leads from terminal block on tube cable bracket. Record and tag lead connections.
5. Remove cable clamp securing control cable to cable bracket by removing three socket head cap screws.
6. Support the X-ray tube and remove the three hex socket flat head screws securing each end of the tube to mounting plates on the yoke assembly.
7. Remove X-ray tube assembly from yoke assembly

2-2 X-Ray Tube Installation

1. Support X-ray tube assembly in a position between mounting plates on yoke assembly. See Illustration 2-1.
2. Attach each end of tube to mounting plates with three hex socket flat head screws. Torque to 50 pound-inches (5.6 N-m).
3. Position control cable in tube cable bracket. Replace cable clamp and secure in place using three socket head cap screws.
4. Connect control cable leads to terminal block on tube cable bracket in accordance with connection data recorded and tagged on leads.

Note: If installing a new tube, use color code and match leads.

5. Position back cover on X-ray tube assembly and secure to tube with four 10-32 $\times$ 1/2 inch button head screws.
6. Cover anode and cathode HV connectors with Silicone Insulating Compound, 46-125224P3. Refer to Direction 46-013288 *Bleeder, High Voltage, Dual Type*.
7. Insert anode and cathode connectors into X-ray tube and secure firmly in place.
8. Check that tube rotates smoothly in yoke assembly when rotated forward and backward between detent positions on yoke mounting plates.
9. Install collimator on X-ray tube. Refer to Section 1-3, steps 1 through 9.
10. If a new tube is installed, the field light must be checked for proper positioning. Refer to Section 1-6, *Adjusting the Field Light*, for this procedure.

Note: If the tube being installed is a replacement, perform "Beam Quality (Half Value Layer)" in Direction 2173227-100 *AMX-4+ Periodic Maintenance*.

SECTION 3 – SERVICING THE HORIZONTAL ARM

This section provides procedures covering the removal and installation of the field serviceable assemblies and parts comprising the horizontal arm assembly.

3-1 Yoke Assembly Removal

CAUTION

Secure the horizontal arm before removing components from the arm. Removing weight allows the counterweight to descend rapidly and the arm to ascend rapidly possibly causing damage to the unit and personal injury. Use Counterweight Lift Tool 46-303491G1.


1. Remove the collimator. Refer to Section 1.
2. Remove the X-ray tube. Refer to Section 2.

WARNING

POTENTIAL FOR EYE INJURY: WEAR EYE PROTECTION WHEN REMOVING THE YOKE ASSEMBLY, SOME PARTS ARE SPRING LOADED AND MAY CAUSE INJURY.

3. Holding cover in place, remove four capscrews securing cover to end of horizontal arm shaft. A spring loaded detent positioning assembly is retained by the cover. Reach behind the cover, hold the detent assembly in place and carefully remove cover. See Illustration 3-1.

ILLUSTRATION 3-1
YOKE ASSEMBLY



4. Carefully remove the detent assembly consisting of two spring covers, two springs and one shaft with one bearing from recessed location in yoke mounting plate.
5. Slide three spring washers, disk, thrust washer and yoke mounting plate from the shaft on horizontal arm.
6. Remove thrust bearing and thrust washer from horizontal arm shaft.

3-2 Yoke Assembly Installation



EYE INJURY! WEAR EYE PROTECTION WHEN REMOVING THE YOKE ASSEMBLY, SOME PARTS ARE SPRING LOADED AND MAY CAUSE INJURY.

1. Lubricate thrust bearing with small amount of lithium (Fiske 630-A-A or equal) grease.
2. Install thrust washer and thrust bearing on horizontal arm shaft. See Illustration 3-1.
3. Place yoke mounting plate with tube supports on horizontal arm shaft.
4. Slide thrust washer, disk, and three spring washers onto horizontal arm. (Concave side of spring washer facing you)
5. Place two springs and spring covers into recessed position in face of mounting plate.
6. Place shaft with bearing on the two springs and carefully install cover over end of horizontal arm shaft.
7. Hold cover in place on horizontal arm shaft to prevent detent assembly from springing loose and secure with four capscrews.
8. Install the X-ray tube. Refer to Section 2.
9. Install the collimator. Refer to Section 1.


3-3 Yoke Mounting Plate Assembly Removal

1. Remove collimator. Refer to Section 1.
2. Remove the X-ray tube. Refer to Section 2.

Right Mounting Plate Assembly


3. Remove decal from cover on end of right tube support to access the cover mounting screws.
4. Remove four screws holding cover to mounting plate shaft.
5. Remove cover, two spring washers, two thrust washers and thrust bearing from mounting plate shaft.
6. Slide shaft, disk, and mounting plate from right tube support.
7. Remove two small bearings and their shaft from right tube support.
8. Remove four screws holding shaft and disk to mounting plate.
9. Remove shaft and disk from mounting plate.

ILLUSTRATION 3-2
RIGHT MOUNTING PLATE ASSEMBLY



Left Mounting Plate Assembly

10. Remove decal from cover on end of left tube support to access the cover mounting screws.
11. Remove four screws holding cover to mounting plate shaft.
12. Remove two spring washers, disk and thrust washer from mounting plate shaft.
13. Slide shaft, thrust washer, dowel pin, disk, and mounting plate from left tube support.
14. Remove thrust washer and dowel pin from shaft.
15. Remove four screws holding shaft and disk to mounting plate.
16. Remove shaft and disk from mounting plate.
17. Remove eight hex socket flat head capscrews holding right and left tube support arms to mounting plate.
18. Remove tube supports.

ILLUSTRATION 3-3
LEFT MOUNTING PLATE ASSEMBLY

3-4 Yoke Mounting Plate Assembly Installation

Right Mounting Plate Assembly

1. See Illustration 3-2. Place mounting plate shaft and disk on mounting plate.
2. Secure shaft and disk to mounting plate with four capscrews.
3. Replace two small bearings and their shaft in right tube support.
4. Slide mounting plate shaft with disk and mounting plate into right tube support.
5. Replace two thrust washers and thrust bearing, two spring washers and cover on mounting plate shaft. (Spring washers, concave side facing out)
6. Secure cover to mounting plate shaft with four capscrews.
7. Replace decal on cover.

Left Mounting Plate

8. See Illustration 3-3. Place mounting plate shaft and disk on mounting plate.
9. Secure shaft and disc to mounting plate with four capscrews.
10. Replace washer and dowel pin on mounting plate shaft.
11. Slide mounting plate shaft with washer, dowel pin, disk and mounting plate into left tube support.

12. Replace washer, disk and two spring washers on mounting plate shaft. (Spring washers, concave side facing out)
13. Secure cover to mounting plate shaft with four capscrews.
14. Replace decal on cover.
15. Install the X-ray tube. Refer to Section 2.
16. Install the collimator. Refer to Section 1.

3-5 Horizontal Arm Assembly Removal


CAUTION

Equipment Damage! Secure the carriage before removing components from the arm. Removing weight allows the arm to ascend rapidly possibly causing damage to the unit and personal injury.


1. Check that the carriage is securely held in place, and that the horizontal arm cannot move upward.
2. Support the X-ray tube and collimator using chain hoist or other suitable equipment.
3. Remove the eight hex socket flat head screws that hold the yoke with tube head and collimator to the mounting plate on the horizontal arm. See Illustration 3-4.

ILLUSTRATION 3-4

REMOVAL AND REMOUNTING OF X-RAY TUBE AND COLLIMATOR ASSEMBLY



4. Lower tube and collimator gently to the floor.
5. Remove seven screws securing horizontal arm assembly to carriage in vertical column. See Illustration 3-5. Discard screws.
6. Remove horizontal arm from carriage.

ILLUSTRATION 3-5
HORIZONTAL ARM MOUNTING

3-6 Horizontal Arm Assembly Installation

Note: The counterweight must be at the bottom of the column or otherwise immobilized.

1. Mount horizontal arm assembly on carriage assembly flanges. See Illustration 3-5.
2. While supporting the arm, secure it firmly in place with six mounting screws. Torque to 110 pound-inches (12.3 N-m).
3. Route brake cable along the anode cable to strap point on side of the column. Secure brake cable to the anode cable with ty-raps.
4. Drape the cable from the side of the column to the terminal box on the back of the column. Insert the brake cable with a Heyco bushing up through the bottom of the terminal box.
5. Connect brake cable leads to terminal block in accordance with connection data recorded and tagged on leads.
6. Position back cover on terminal box and secure with two screws.

CAUTION


The horizontal arm must be secured before installing components on the arm. Do not rely on latch alone to hold it in place. With weight removed from the arm, the counterweight can descend rapidly and the arm ascend rapidly possibly causing damage to the unit and personal injury. Use Counterweight Lift Tool 46-303491G1.

7. Firmly grasp the horizontal arm and using your weight to counterbalance the counterweight, approximately 75 pounds (34 kg), carefully move it down to the top of the unit and into parked position. Do not rely on latch alone to hold the arm in place. Check that the arm is locked down before slowly removing your weight and releasing it.
8. Re-mount the X-ray tube and collimator assembly per Illustration 3-4.
9. Check that the arm moves freely up and down the column.

3-7 Horizontal Arm Brake and Cable Assembly Removal

1. Remove horizontal arm assembly. Refer to Section 3-5.
2. Position horizontal arm assembly with latch block up.
3. Extend Horizontal Arm until clearance hole on underside is exposed.
4. Remove cable clamp. See Illustration 3-6.
5. Disconnect Brake Cable connection.
6. Remove the two stop block capscrews securing the brake and cable assembly to the end of the outer tube. See Illustration 3-6.

ILLUSTRATION 3-6
HORIZONTAL ARM OUTER, CENTER AND INNER TUBES



7. Return arm to upright position.
8. Slide the center and inner tubes from the outer tube.
9. Place center and inner tubes in upside down position.
10. Remove four screws holding latch block to center tube. Remove block. This provides access to the mounting plate holding screw. See Illustration 3-6.
11. Remove mounting plate capscrew holding brake and cable assembly to the center tube.
12. Remove the two brake capscrews holding the brake & cable assembly to the rear of the center tube.
13. Remove two or three screws (depends on unit you have) holding the connector bracket and ground wire to the center tube.
14. Remove two retaining rings from inner shaft on center tube. Slide shaft out permitting removal of two spacers and two bearings.
15. Carefully slide the brake and cable assembly, and the inner tube from the center tube.
16. Remove two capscrews holding the brake and cable assembly bracket to the inner tube.
17. Remove brake and cable assembly from inner tube.

3-8 Horizontal Arm Brake and Cable Assembly Installation


1. Mount the brake and cable assembly on the inner tube.
2. Secure the brake and cable assembly bracket to the inner tube with two capscrews. See Illustration 3-6.
3. Carefully slide the brake and cable assembly, and the inner tube into the center tube.
4. Replace the inner shaft, two spacers and two bearings. Secure with two retaining rings.

5. Place center and inner tubes in an upside down position.
6. Replace the two brake capscrews securing the brake & cable assembly to the rear of the center tube.
7. Secure the connector bracket and ground wire to the center tube with two or three screws (depends on unit you have).
8. Replace mounting plate capscrew holding brake and cable assembly to the center tube.
9. Replace the latch block and secure it with four screws to center tube.
10. Place the outer tube in an upside down position.
11. Slide the center and inner tubes into the outer tube.
12. Replace the two stop block capscrews and secure the brake and cable assembly to the end of the outer tube.
13. Connect Brake Cable.
14. Extend Horizontal Arm until clearance hole on underside is exposed.
15. Replace cable clamp.
16. Return arm to upright position.
17. Replace horizontal arm assembly. Refer to Section 3-6.

3-9 Horizontal Arm Brake Cable Removal

1. Remove horizontal arm brake and cable assembly. Refer to Section 3-7.
2. Remove cable clip from bar support assembly by removing one machine screw. See Illustration 3-7.

ILLUSTRATION 3-7
BRAKE AND CABLE ASSEMBLY



3. Remove anchor nut from brake cable by turning counterclockwise. This nut has been secured with set screws and Loctite.
4. Remove cable anchor end from bar support.
5. Slide ball end of cable from slot in cable bracket and remove from path around pulley end on mounting plate.
6. Remove the brake cable from brake and cable assembly.

3-10 Horizontal Arm Brake Cable Installation

1. Slide ball end of brake cable into slot in cable bracket on end of chain on brake and cable assembly. See Illustration 3-7.
2. Route anchor end of cable around pulley on end of mounting plate to opposite end of support bar and insert in slot.

3. Replace anchor nut on cable end and take up slack by turning nut clockwise. Do not overtighten.
4. Replace cable clip, with brake cable in place, on bar support assembly and secure with one machine screw.
5. Check that the cable and chain move freely.
6. Loctite the anchor nut.
7. Install horizontal arm brake and cable assembly. Refer to Section 3-8.

3-11 Horizontal Arm Brake Removal

1. Remove horizontal arm brake and cable assembly. Refer to Section 3-7.
2. Disconnect brake electrical cable connectors to remove from connector bracket.
3. Loosen remaining capscrew holding brake to mounting plate. See Illustration 3-7.
4. Slide brake off brake sprocket shaft.

3-12 Horizontal Arm Brake Installation

1. Slide brake onto brake sprocket shaft.
2. Secure brake with single capscrew in end of mounting plate. See Illustration 3-7.
3. Position brake electrical cable connectors in connector bracket and reconnect.
4. Install horizontal arm brake and cable assembly. Refer to Section 3-8.


SECTION 4 – SERVICING THE VERTICAL COLUMN

This section provides procedures covering the removal and installation of the field serviceable assemblies subassemblies and parts comprising the vertical column assembly.

4-1 Remove Vertical Column Assembly

1. Remove x-ray tube, collimator, and horizontal arm assemblies per Section 3-5.
2. Disconnect control cable leads from terminal block. Record and tag leads. See Illustration 4-1.
3. Remove Heyco bushing and cable through bottom of terminal box.
4. Clear all cables and wires from the column area.
5. Remove vertical column cap by removing two hex socket button head capscrews. See Illustration 4-1.
6. Attach lifting tool (46-303362P1) to top of column with 2 screws.
7. Attach chain hoist to lifting tool (46-303362P1) and take up slack. The chain hoist must be capable of lifting a seven foot (2.1 meter), 175 lb. (80 kg) column a minimum of two feet (61 cm), preferably 4 feet (122 cm)..


ILLUSTRATION 4-1
COLUMN TRIM COVER REMOVAL



8. Remove trim cover from bottom of vertical column by removing three screws, two clip nuts and two binding head screws. See Illustration 4-1.
9. Remove 5 bearing plate screws and washers securing column to the column support base.
10. Remove two screws securing brake cable cover and disconnect the lower brake connectors. See Illustration 4-2.

11. Move brake cable connectors and cover away from column support base.
12. Break loose the column support bearing from the column support base by jacking the column out using three 1/4-20 screws thru the base plate and against the support base.
13. Lift the column and bearing plate out of the base.
14. Lower the column onto the floor or a table capable of supporting its weight. Carriage side must be upwards.
15. Remove the chain hoist and lifting tool.

ILLUSTRATION 4-2
VERTICAL COLUMN REMOVAL



4-2 Install Vertical Column Assembly

1. Attach lifting tool (46-303362P1) to top of pulley and brake assembly with two capscrews. See Illustration 4-2.
2. Attach chain hoist to lifting tool.

Note: Make sure screw holes are clean of oil, grease or debris.

3. Lower the column into the base. Carefully seat the bearing plate in the base with the mounting holes lined up and the lower brake connectors opposite the opening in the base. See Illustration 4-2.
4. Apply Loctite 242 (medium strength) to threads of five new column bearing plate capscrews.

Note: Torque cap screws to 100 pound-inches (11.3 N-m).

5. Install five new capscrews and five new flat washers and secure bearing plate to the column support base. See Illustration 4-2.
6. Connect brake connectors to brake. See Illustration 4-2.
7. Secure cover on column support base with two screws.
8. Replace trim cover at base of column. Refer to Section 5.
9. Insert control cable with Heyco bushing through bottom of terminal box located on back of vertical column. See Illustration 4-1.
10. Connect control cable leads to terminal block as recorded on tags.
11. Remove the chain hoist and lifting tool.
12. Replace trim cover on top of vertical column and secure with two hex socket button head capscrews. Refer to Section 5.
13. Replace horizontal arm. Refer to Section 3-6.
14. Then continue in Section 3-6 to remount x-ray tube and collimator assembly.
15. Check counterweight balancing. Refer to Section 10.

4-3 Counterweight Cables

Note: The follower cable should be replaced only if it shows signs of wear or deterioration, and can only be removed using the procedure in Sections 4-3-5 and 4-3-6.

The main and safety cable should be replaced periodically as defined in the Planned Maintenance section of the Operator Manual, or where elsewhere directed.

CAUTION


This procedure requires that the FE work at the top of the vertical column. Ensure proper safety precautions are taken against falls or loss of balance.

This procedure can result in dropping the counterweight to the bottom of the column if done incorrectly. If this occurs use the procedure in Sections 4-3-5 and 4-3-6 to replace the cables and reattach the counterweight.

1. Set the main power circuit breaker to OFF.
2. Position horizontal arm and tube at least 90 degrees out of park, so that it can be positioned at lowest possible position.
3. Raise horizontal arm to a comfortable working height.

Note: Use a C-clamp to hold horizontal arm in position during this procedure if desired. Protect the column from clamping damage if this is done.

4. Remove column cap.
5. Remove four binding head screws holding cover plate on the carriage assembly. See Illustration 4-3.

ILLUSTRATION 4-3
CABLE REPLACEMENT


4-3-1 Main Cable

1. Loosen locking nut on top of stud end of main cable. Main cable is in front of follower cable.

CAUTION

Sudden and unexpected motion can occur: After next step the counterweight is held only by the safety cable.

2. Remove holding nut on bottom of stud end of main cable.
3. Remove cable stud from carriage cable mounting strut.
4. Move horizontal arm to bottom of column. This raises the counterweight to the top of the column.
5. Remove counterweight trim weights (if present). See Illustration 4-4.



6. Remove screws holding cable retainer clip from top of counterweight using a magnetic nut driver or screwdriver.

Note: The positioning is awkward and the screw can be easily dropped.

7. Push main cable ball end out from counterweight keyhole slot.
8. Pull main cable stud end through top of column.
9. Remove old cable.
10. Push new main cable stud end through top of column over top of upper pulley.
11. Push ball end into counterweight slot.
12. Install cable retainer clip.
13. Thread new main cable through back of carriage.
14. Raise vertical arm to comfortable working height.
15. Apply Loctite to threads of cable stud. Install jam nut on cable stud.
16. Install cable stud through carriage horizontal bar and install holding nut.
17. Set cable to tension in accordance with Section 4-3-4. Tighten cable jam nut.

Note: Do not allow cables to twist while setting tension.

4-3-2 Safety Cable

1. Loosen locking nut on top of stud end of safety cable. Safety cable is not in front of follower cable.



After next step the counterweight safety mechanism is temporarily inactive until the new safety cable is installed.

2. Remove holding nut on bottom of stud end of safety cable.
3. Remove cable stud from carriage cable mounting strut.
4. Move horizontal arm to bottom of column. This raises the counterweight to the top of the column.
5. Remove screws holding cable retainer clip from top of counterweight using a magnetic nut driver or screwdriver.

Note: The positioning is awkward and the screw can be easily dropped.

6. Push safety cable ball end out from counterweight groove.
7. Pull safety cable stud end through top of column.
8. Remove old safety cable from column.
9. Push new safety cable stud end through top of column over top of upper pulley.
10. Push ball end into counterweight slot.
11. Install cable retainer clip.
12. Install trim weights.
13. Thread new safety cable through back of carriage.
14. Raise vertical arm to comfortable working height.
15. Apply Loctite to threads of cable stud. Install jam nut on cable stud.
16. Install cable stud through carriage horizontal strut and install holding nut.
17. Set cable tension in accordance with Section 4-3-4. Tighten cable jam nut.

Note: Do not allow cables to twist while setting tension.

4-3-3 Reassemble Unit

1. Replace carriage access plate and four binding head screws.
2. Install column cap.
3. Move horizontal arm through full range of motion, ensuring that it moves smoothly with equal force throughout its range. Ensure that both cables run smoothly through upper pulley.

4-3-4 Vertical Column Cable Tension

Tension in the vertical column cables can be measured two ways: when the column is horizontal (as when replacing the counterweight) or vertical (as when replacing the main and safety cables). Table 4-1 lists the tension values required for the different attitudes.

TABLE 4-1
VERTICAL COLUMN TENSION


COLUMN ATTITUDE	CABLE TENSION		
	MAIN	SAFETY	FOLLOWER
HORIZONTAL	60 LB (27 KG)	5-10 LB (2.3-4.5 KG)	N/A
VERTICAL	N/A	N/A	10-15 LB (4.5-6.8 KG)
N/A – Not applicable; value not measured			

4-3-5 Remove Follower Cable

If replacement of follower cable is required, all three cables should be replaced to ensure that the periodic maintenance time interval specified for replacement of the cables remains the same.

Note: Use this procedure for follower cable service, or in the event the procedure in Section 4-3 cannot be successfully performed.

1. Remove vertical column from unit. Refer to Section 4-1. Vertical Column removal is necessary in order to replace a follower cable.
2. Remove cover plate from the carriage assembly by removing four binding head screws holding it in place. See Illustration 4-3.

ILLUSTRATION 4-5
CABLE REPLACEMENT

3. Move carriage near to bottom of column. This also brings counterweight to top of column where it can be removed.
4. Loosen locking nuts on top of stud ends of main and safety cables.
5. Remove holding nuts on bottom of stud ends of main and safety cables.
6. Remove cable studs from carriage cable mounting strut.
7. Push ball end of follower cable out of slot in carriage cable mounting strut.
8. Move carriage toward top of column until ball end of follower cable is exposed. Move ball end of cable to bottom of column.
9. Remove follower cable from lower pulley and leave inside bottom of column.
10. Remove pulley and brake assembly. Refer to Section 4-5.
11. Remove counterweight from top end of column. Refer to Section 4-7.

Note: At this point all three cables are detached from the counterweight and can be removed from the vertical column assembly.

12. Remove main, safety and follower cables.

4-3-6 Install Follower Cable

1. Install counterweight assembly. Refer to Section 4-8.
2. Install upper pulley and brake assembly. Refer to Section 4-6.
3. Reach in the bottom of the column, retrieve the follower cable and position it around the lower pulley. See Illustration 4-5.
4. Place the ball end of the follower cable in the small square cut-out in the column.
5. Align the carriage cable mounting strut over the base end of the follower cable.
6. Place ball end of follower cable in slot on carriage cable mounting strut.
7. Replace top locking nuts on stud ends of main and safety cables.
8. Insert main safety cable studs into carriage cable mounting strut holes.
9. Replace holding nuts on bottom of stud ends of main and safety cables. Do not tighten at this time.
10. Make certain all the cables are still in their pulley grooves. (Follower cable is on right side as you look at it.)
11. Tighten the main cable until there is no slack in the follower cable.
12. Set follower cable to tension value found in Section 4-3-4. Tighten main cable jam nut.
13. Apply Loctite and tighten the locking nuts.
14. Check that the carriage moves freely up and down the column.
15. Replace cover plate on the carriage assembly and secure with four binding head screws. See Illustration 4-2.
16. Move carriage to top of column. This moves counterweight to bottom of column.
17. Install vertical column. Refer to Section 4-2.

4-4 Upper Brake, Pulley, and Bearings


4-4-1 Disconnect Brake

1. Remove the vertical column. Refer to Section 4-1.
2. Tag either the main or safety cable. This will eliminate confusion when securing them in their correct positions on the carriage mounting strut.
3. Detach the main and safety cables from the carriage cable mounting strut. Refer to Section 4-3.
4. Remove two screws in the bottom of terminal box located on back of vertical column and lift off cover.
5. Disconnect upper brake cable leads from terminal block. (It is important to tag and record terminal points, otherwise there is a possibility of blowing internal diode in brake if leads are improperly installed.)
6. Remove brake cable from terminal box and channel on back of column.

4-4-2 Remove Upper Pulley and Brake

7. Remove capscrews securing pulley and brake assembly to top of column. See Illustration 4-6.
8. Remove Heyco bushing from column where brake cable enters pulley and brake assembly. Bushing extends into pulley and brake casting.


ILLUSTRATION 4-6
UPPER PULLEY AND
BRAKE ASSEMBLY



9. Slowly withdraw pulley and brake assembly from top of column while sliding the brake cable through slot in top of column. These must be done together or the brake cable will be severed.
10. Pull main and safety cables from the column.
11. Remove the nuts from the stud ends of both cables.
12. Remove each cable from the pulley and brake assembly by first lifting it from its pulley groove, moving it to the side of the pulley and then pulling it free. The ball ends of these cables still remain attached to counterweight.
13. Remove the pulley and brake assembly.

4-4-3 Replace Upper Pulley and Bearings

1. Remove upper pulley and brake. Refer to Sections 4-4-1 and 4-4-2.
2. Place a tag on either the main cable or safety cable. This will aid in re-installing them later in their correct positions.
3. Remove pulley set screws. See Illustration 4-7.

ILLUSTRATION 4-7
UPPER BRAKE AND PULLEY

4. Remove retaining ring from end of pulley shaft.
5. Remove pulley shaft bearing, key, pulley and cables from mounting bracket.
6. Remove bearing from mounting bracket.
7. Place bearing on pulley shaft and insert key.
8. Place pulley in position in the mounting bracket.
9. Insert pulley shaft, bearing and key into mounting bracket and through pulley.
10. Place bearing in mounting bracket on end of pulley shaft.
11. Place retaining ring on end of pulley shaft.
12. Replace pulley set screws and secure firmly.

4-4-4 Install Upper Pulley, Bearings and Brake

1. Place pulley and brake assembly near upper end of vertical column. See Illustration 4-6.
2. Remove the nuts from the stud ends of main and safety cables.
3. Insert each cable into pulley and brake assembly along the side of the pulley then lift it into its assigned pulley groove. The ball ends of these cables are attached to the counterweight.
4. Insert main and safety cables into the column through slot provided.
5. Check that the brake cable is exiting pulley and brake assembly through hole in casting.
6. Replace pulley and brake assembly in top of column while sliding the brake cable through slot in top of column. This must be done together or the brake cable may be damaged.
7. Replace the Heyco bushing in column where brake cable passes through column and enters pulley and brake assembly. Bushing extends into pulley and brake casting.

Note: Torque cap screws to 50–70 pound–inches (5.6 to 7.8 N–m).

8. Apply Loctite and replace cap screws securing pulley and brake assembly to top of column.

4-4-5 Connect Brake

1. Insert brake cable through top of channel on back of column and down into terminal box.
2. Connect upper brake cable leads to terminal block. See Illustration 4–2.
3. Replace terminal box cover located on back of vertical column and secure with two screws.
4. Install main and safety cables on the carriage cable mounting strut. Refer to Section 4–4.
5. Install vertical column. Refer to Section 4–2.

4-4-6 Remove Upper Pulley, Bearings and Brake (Alternate)

Note: This procedure can be used to replace the brake without removing the vertical column

1. Remove horizontal arm. Refer to Section 3–5. It is not necessary to remove yoke from arm.
2. Disconnect upper brake. Refer to Section 4-4-1.
3. Disconnect main and safety cables from carriage. Refer to Section 4-6.
4. Move carriage to bottom of column.
5. Remove four capscrews securing pulley and brake assembly to top of column.
6. Remove the Heyco bushing in column where brake cable passes through column and enters pulley and brake assembly. Bushing extends into pulley and brake casting.
7. Remove pulley and brake assembly from top of column while sliding the brake cable through slot in top of column. This must be done together or the brake cable may be damaged.
8. Remove screws holding upper brake to the mounting bracket.


4-4-7 Install Upper Pulley, Bearings and Brake (Alternate)

Note: This procedure can be used in conjunction with Section 4-4-6.

1. Insert screws holding upper brake to the pulley and brake assembly casting and secure in place.
2. Check that the brake cable is exiting pulley and brake assembly through hole in casting.
3. Replace pulley and brake assembly in top of column while sliding the brake cable through slot in top of column. This must be done together or the brake cable may be damaged.
4. Replace the Heyco bushing in column where brake cable passes through column and enters pulley and brake assembly. Bushing extends into pulley and brake casting.

Note: Torque cap screws to 50–70 pound–inches (5.6 to 7.8 N–m).

5. Apply Loctite and replace capscrews securing pulley and brake assembly to top of column.
6. Connect brake. Refer to Section 4-4-5.
7. Move carriage to top of column and insert main and safety cables in carriage cable mounting strut. See Illustration 4–8.

ILLUSTRATION 4-8
CARRIAGE ASSEMBLY

8. Tighten the main and safety cable stud end nuts and lock in place with the holding nuts.
9. Set cable tension in accordance with Section 4-3-4. Tighten cable jam nuts.
10. Replace cover plate on the carriage assembly and secure with four binding head screws. See Illustration 4-3.
11. Replace horizontal arm with yoke. Refer to Section 3-6.

4-5 Counterweight


4-5-1 Remove Counterweight

1. Remove vertical column. Refer to Section 4-1.
2. Remove the main, safety and follower cables. Refer to Section 4-3.
3. Remove upper pulley and brake assembly. Refer to Section 4-4.



PERSONAL INJURY! THE COUNTERWEIGHT WEIGHS APPROXIMATELY 100 POUNDS (45 KG). USE THE PROPER SAFETY PRECAUTIONS WHEN HANDLING THE COUNTERWEIGHTS.

4. Slide counterweight from vertical column.
5. Remove top cable retainer clip from counterweight assembly by removing two slotted hex washer head screws. See Illustration 4-9.

ILLUSTRATION 4-9
COUNTERWEIGHT

6. Remove main and safety cables from counterweight assembly by dislodging ball ends of cables from retaining notch in counterweight.
7. Remove bottom follower cable retainer clip from counterweight assembly by removing two slotted hex washer head capscrews.
8. Remove follower cable from counterweight assembly by dislodging ball end of cable from retaining notch in counterweight.

4-5-2 Install Counterweight

1. Replace main and safety cables in counterweight assembly by inserting ball ends of cables into retaining notch in top of counterweight. See Illustration 4-9.
2. Replace cable retainer clip in top of counterweight assembly and secure with two slotted hex washer capscrews.
3. Place follower cable ball end into retaining notch in bottom of counterweight. The follower cable must be installed on the right side to function properly.
4. Replace follower cable retainer clip in bottom of counterweight assembly and secure with two slotted hex washer head capscrews.
5. Slide follower cable down inside of vertical column as far as possible.
6. Slide counterweight assembly into vertical column.
7. Check counterweight balancing. Refer to Section 10.

4-5-3 Counterweight Balancing

Note: The column assembly trim cover must be removed from top of column assembly as shown in Section 4-1.


1. Secure the collimator hand switch closed with a ty-wrap or piece of tape to energize the vertical column brake (as well as all other brakes and latch solenoid).
2. With all hardware in place, check the force required to move the tube/arm assembly upward and downward. If the force needed to move it up is equal to the force needed to move it down, the vertical column is balanced. Proceed to Step 5.
3. If the vertical moving forces are unequal, eight oz. (0.23 kg) counterweights (part number 46-279784P1) should be added or removed (approximately four are installed at the factory).
 - a. Move the tube/arm to its lowest possible position. This will raise the counterpoise mass and counterweight to the top of the vertical column where the counterweight can be accessed.
 - b. Loosen and remove the two locknuts and washers that hold the counterweights to the counterpoise.
 - c. Add or remove counterweights as required,
 - d. Replace the washers and locknuts, and tighten locknuts.
4. Repeat Step 2. If the moving forces are equal, proceed to Step 5. If they are still unequal, repeat Step 3.
5. Replace trim cover on top of vertical column and secure with two hex socket button head capscrews.

4-6 Carriage Assembly

4-6-1 Remove Carriage

1. Remove upper pulley and brake assembly. Refer to Section 4-5.
2. Remove follower cable end from carriage strut by loosening its holding nut. See Illustration 4-10.
3. Slide carriage to top of column and remove it.

ILLUSTRATION 4-10
CARRIAGE ASSEMBLY



4-6-2 Install Carriage


1. Slide carriage assembly into top of column. Check that horizontal arm mounting plate is at top of carriage.
2. Install pulley and brake assembly. Refer to Section 4-6.
3. Place follower cable end into carriage strut and secure in place by tightening its holding nut. See Illustration 4-10.

4-6-3 Remove Carriage Bearings

1. Remove carriage assembly. Refer to Section 4-6-1.

2. Remove four hex head shoulder bolts securing four forward facing bearings to side flanges on carriage. See Illustration 4-11.
3. Remove carriage side facing bearing shoulder bolts and bearings from carriage.

ILLUSTRATION 4-11
CARRIAGE BEARINGS




4-6-4 Install Carriage Bearings

1. Lubricate bearings if necessary.
2. Attach four forward facing bearings to carriage mounting flanges with hex head shoulder bolts. See Illustration 4-11. Apply Loctite and torque to 40-60 pound-inches (4.5 to 6.7 N-m)
3. Install side facing bearings at mounting locations on each side of carriage with shoulder bolts.
4. Install carriage assembly in vertical column assembly. Refer to Section 4-14.

4-7 Lower Brake Removal

1. Remove vertical column. Refer to Section 4-1.
2. Remove retaining ring from bottom of column base shaft. See Illustration 4-12.

ILLUSTRATION 4-12
LOWER BRAKE AND BEARINGS



3. Pull lower column bearing from column base shaft.
4. Remove screws and washers securing lower brake to column bearing plate.
5. Remove lower brake from column base shaft.

4-8 Lower Bearings and Baseplate Removal

1. Remove vertical column. Refer to Section 4-1.
2. Remove lower brake. Refer to Section 4-7.
3. Remove sleeve and key from column base shaft.
4. Remove upper column bearing and column bearing plate from column base shaft.

Note: It is not normally necessary to take the column bearing plate apart, but the stop disk can come out of its position in the detent ring.

4-9 Lower Brake and Bearings Installation


Note: Install stop disk in column bearing plate if it has come out.

1. Install upper column bearing in column bearing plate. See Illustration 4-12.
2. Place column bearing plate on column base shaft. Ensure that column bearing plate sits flat on bottom of column.
3. Install sleeve and key on column base shaft.
4. Place lower brake on column base shaft.
5. Secure lower brake to bearing plate with screws and washers. Torque to 40–60 pound-inches (4.5 to 6.7 N-m).
6. Slide lower bearing on column base shaft and secure with retaining ring.
7. Install vertical column assembly. Refer to Section 4-2.

4-10 Column Base Removal

1. Remove vertical column assembly. Refer to Section 4-1.
2. Remove follower cable. Refer to Section 4-3-5.
3. Remove lower brake and bearings. Refer to Sections 4-7 and 4-8.
4. Remove hex head bolt and flat washers securing column base assembly to vertical column. See Illustration 4-13.

ILLUSTRATION 4-13
COLUMN BASE



5. Pull column base assembly and spacers from vertical column.

4-11 Column Base Installation


1. Insert column base assembly into vertical column. Line up holes with those in vertical column. See Illustration 4-13.
2. Position spacers so that holes line up with column base and column.

3. Apply Loctite to the hex head bolts. Secure column base assembly to vertical column with bolts and flat washers. Torque bolts to 30–35 pound-feet (41 to 47 N-m).
4. Install lower brake and bearings. Refer to Sections 4-7 and 4-8.
5. Install follower cable. Refer to Section 4-3-5.
6. Install vertical column assembly. Refer to Section 4-2.

4-12 Lower Pulley and Bearings Removal

1. Remove column base. Refer to Section 4-10.
2. Loosen two set screws on pulley hub. See Illustration 4-14.

ILLUSTRATION 4-14
LOWER PULLEY AND BEARINGS




3. Remove external retaining rings from pulley ends of pulley shaft.
4. Pull out shaft and remove pulley.
5. Remove bearings from column base casting.

4-13 Lower Pulley and Bearings Installation

1. Insert bearings in bearing seats in column base casting. See Illustration 4-14.
2. Push pulley shaft through left bearing. Slide pulley on shaft with pulley hub toward right bearing.
3. Insert pulley shaft in right bearing and secure shaft with one external retaining ring in groove on each end of shaft.
4. Tighten two set screws on pulley hub.
5. Install column base. Refer to Section 4-2.

4-14 Column Support Assembly Removal

1. Remove vertical column assembly. Refer to Section 4-1.
2. From beneath the bare assembly, remove capscrews and washers securing column support to bottom of base assembly. See Illustration 4-15.

ILLUSTRATION 4-15
COLUMN SUPPORT ASSEMBLY

3. Remove hex head capscrews, lockwashers and flat washers securing the column support to the base assembly.
4. Remove column support from base assembly.

4-15 Column Support Assembly Installation

1. Place column support in base assembly. Orient opening in side of column support so it is opposite terminal board in base assembly. Line up mounting holes. See Illustration 4-15.
2. Secure the column support to the base assembly with hex head capscrews, lockwashers and flat washers. Torque to 100–120 pound-inches (11.2 to 13.4 N-m).
3. Secure bottom of column support to bottom of front tongue in base assembly with cap screws and washers. Torque to 100–120 pound-inches (11.2 to 13.4 N-m).
4. Install vertical column assembly. Refer to Section 4-2.


SECTION 5 – SERVICING THE CASSETTE DRAWER, COVERS, HORIZONTAL ARM LATCH AND HANDSWITCH

This section provides procedures covering the removal and installation of trim covers.

5-1 Cassette Drawer Removal

1. Remove rear cover from unit by removing two screws securing it to base assembly. See Illustration 5-1.
2. Remove five binding head screws and flat washers securing cassette drawer hinge to base assembly.
3. Pull cassette drawer to open position. It is retained by magnetic catch. Lift drawer up and out of base assembly.

ILLUSTRATION 5-1
CASSETTE DRAWER



5-2 Cassette Drawer Installation


1. Place cassette drawer in position on base assembly.
2. Attach cassette drawer hinge to base with five binding head screws and flat washers.
3. Install rear cover to base assembly with two screws.
4. Close cassette drawer. Magnetic catch will hold it in its place.

5-3 Side Covers Removal

Note: This procedure covers the removal of both the right and left trim covers.

1. Turn power off by turning circuit breaker OFF.
2. Remove front cover from unit by pulling corners of cover from ball stud retainer sockets on mainframe. See Illustration 5-2.
3. Remove rear cover by removing two screws retaining it on base assembly.
4. Disconnect hand switch on right side trim cover.
5. Remove two 4-40 x 3/8 inch binding head screws and one 6-32 x 1/2 inch black oxide binding head screw holding each cover to base and mainframe assemblies.
6. Pull covers loose from five ball stud retaining sockets and carefully remove from unit.

ILLUSTRATION 5-2
SIDE COVERS



5-4 Side Cover Installation

1. Turn power off by turning circuit breaker OFF.

Note: This procedure covers the installation of both the right and left trim covers.


2. Install side covers by first pressing the five ball studs into retainer sockets on base and mainframe assemblies. See Illustration 5-2.
3. Secure each side cover to base and mainframe assemblies with two 4-40 x 3/8 inch binding head screws and one 6-32 x 1/2 inch black oxide binding head screw.

4. Connect handswitch on right side trim cover.
5. Install rear cover with two screws.
6. Install front cover on unit by pressing ball studs on corners of cover into retainer sockets.

5-5 Top Cover Removal

1. Move horizontal arm to top of column.
2. Turn power off by turning circuit breaker OFF.
3. Open cassette drawer, reach up into top cover and rotate two spring latches to release cover from the mainframe.
4. Remove three hex socket capscrews holding collar to horizontal arm latch and lift off collar and spacer beneath it. See Illustration 5-3.

ILLUSTRATION 5-3
TOP COVER



5. Lift top cover straight up until right and left hinges lock in vertical position.

Note: With top cover extended on hinges, the cover can be tilted toward the front or rear of the unit to provide access to either area for servicing. If complete removal of the cover is necessary then continue with steps 5 through 9.

6. Disconnect wires from key switch assembly.
7. Disconnect ribbon cable from display at printed wire board.
8. Remove the shouldered screws holding the right and left hinges to the top cover brackets.
9. Remove the top cover.

5-6 Top Cover Installation

1. Move horizontal arm to top of column.
2. Turn power off by turning circuit breaker OFF.
3. Extend right and left hinges from mainframe to vertical position. See Illustration 5-3.


4. Lift top cover and position brackets on hinges.
5. Secure top cover brackets to right and left hinges with two shouldered screws.
6. Connect ribbon cable from display to printed wire board.
7. Connect wires to key switch assembly.
8. Carefully lower cover into position on top of mainframe.
9. Open cassette drawer, reach up into top cover and rotate two spring latches to lock cover to mainframe.
10. Place spacer and collar on horizontal arm latch and secure collar in place with three hex socket capscrews.

5-7 Remove Latch Assembly and Components (Hall Effect Switch)

1. Move horizontal arm out of the way.
2. Shut power off by turning circuit breaker OFF.
3. Open top cover. Refer to Section 5-5.


Note: If you need more access to terminal strip, you may remove four hex nuts attaching low voltage power supply to mainframe and slide power supply backward out of mainframe housing. See Illustration 5-4.

ILLUSTRATION 5-4
LATCH ASSEMBLY (HALL EFFECT LATCH SWITCH)



4. Disconnect leads from terminal strip (see Illustration 5-5):

- > Sensor Board: three leads to quick disconnect connectors. (#1, #2 and #3).
- > Solenoid: two leads to terminal strip. Tag and label terminals. Must reconnect properly on #3 & #4 or will blow diode in solenoid.

ILLUSTRATION 5-5
LATCH WIRING

5. Remove latch assembly from mainframe by removing four hex socket head capscrews.
6. Remove latch block from spacer by removing four hex socket head capscrews.
7. Push spring pin out of latch block.
8. Remove latch.
9. Remove latch spring by removing binding head screw.
10. Remove hex nut from top of solenoid.
11. Remove solenoid from recess in latch block.
12. Remove the two flat head screws in side of latch base.

NOTICE

Use care to avoid electrostatic discharge damage to the Hall effect sensor circuit board.

13. Remove the cover and the Hall effect sensor circuit board.

5-8 Install Latch Assembly and Components (Hall Effect Switch)

1. Shut off power by turning circuit breaker OFF.
2. Install solenoid in recess in latch block. See Illustration 5-4.
3. Install hex nut on top of solenoid.

NOTICE

Use care to avoid electrostatic discharge damage to the Hall effect sensor circuit board.

4. Install Hall effect sensor circuit board and cover with two flat head screws.
5. Attach latch spring to latch with binding head screw.
6. Install latch in latch block.
7. Insert spring pin in latch block to retain latch in position.

8. Install latch block on spacer and fasten with four hex socket head capscrews.
9. Install latch assembly on mainframe and fasten with four hex socket head capscrews.
10. Connect leads from latch assembly:
 - > Hall Effect Sensor Board: three leads to terminal strip #1, #2 and #3.
 - > Solenoid: two leads to terminal strip in accordance with labels on terminals.
11. If you moved low voltage power supply out of the way in Section 5-7, slide power supply onto mainframe housing and fasten with four hex nuts.
12. Close top cover. Refer to Section 5-6.


5-9 Remove And Install Hand Switch Cable

To remove the hand switch cable

1. Disconnect the cable from the connector on the body of the AMX-4+.
2. Remove the 6-32 screw securing the cable clamp, then remove the clamp and star washer.

To install the hand switch cable, perform the above steps in reverse order. Be sure to align the cable to the cover as indicated in Illustration 5-6.


ILLUSTRATION 5-6
HAND SWITCH CABLE



5-10 Install Hand Switch

1. Plug end of coil cord into X-Ray Hand Switch.
2. Place X-Ray Hand Switch on support bracket.

ILLUSTRATION 5-7
X-RAY HAND SWITCH



SECTION 6 – SERVICING THE BATTERIES

⚠ WARNING

PERSONAL INJURY ! REMOVE ALL JEWELRY AND RINGS WHEN HANDLING THE BATTERIES OR REACHING IN THE BATTERY COMPARTMENT. USE EXTREME CAUTION, HIGH AMPERAGE AND VOLTAGE IS PRESENT. USE INSULATED TOOLS FOR WORKING ON BATTERIES. WEAR SAFETY GLASSES.

During removal:


- Reference to the right and left sides of battery compartment are relative to the view shown in Illustration 4.
- Cable numbers (#s) indicated in Illustration 6-3 and in this instruction are marked on the parts.
- Insulator numbers are shown for convenience in description, but are not marked on the parts.
- There are no Insulators #3 and #7.
- Battery numbers may or may not be marked on batteries. Mark batteries with numbers if desired.
- Disconnect cables only in the order instructed, and only when instructed.
- Don't remove paper backing from any insulating separator pad.

6-1 Remove Battery Cover

1. Place the unit on a level floor. Extend the horizontal arm and rotate it to the front end (vertical column end) of the unit. See Illustration 6-1.

Note: Distribution of weight is necessary to remove the battery compartment cover.

ILLUSTRATION 6-1
BATTERY COMPARTMENT ACCESS



2. Verify that AMX key switch is OFF, circuit breaker is OFF, and charger cord is unplugged.
3. Remove battery ground cable marked #2 from the ground stud. Tape to ensure it does not short. SEPARATE IT FROM ALL OTHER CABLES.

Note: All battery lead connectors must be insulated to prevent electrical contact with the frame or other battery terminals.


4. Tape cable marked #5 separately from ALL other cables.

5. Remove all remaining cables from ground stud. See Illustration 6-1. Tape each cable separately. Keep ALL cables separate.
6. Remove battery lead cable marked #4 from terminal on Circuit Breaker in AMX1 A3 Module. Be certain boot covers the connector and it cannot short out. Tape connector to ensure it will not short.
7. Remove battery compartment cover capscrews and cover.
8. Remove insulator board from compartment entrance.

6-2 Disconnect Battery Sense Harness (If Required)

1. If unit has a battery sense harness, disconnect harness from J11 of the battery test terminal strip. See Illustration 6-2.
2. Disconnect sense leads from J9 and J10 on left hand side of the battery test terminal strip.


ILLUSTRATION 6-2
SENSE LEADS ON BATTERY TEST TERMINAL STRIP (46-321370G1)
PRINTED WIRING BD ASM
46-321370G1



6-3 Disconnect and Remove Batteries

1. Disconnect Cable #2 from Battery #9 (-) and slide boot over cable terminal. See Illustration 6-3.
2. Disconnect Cable #5 from Battery #9 (+), and slide boot over cable terminal.
3. Remove Battery #9, Insulator #8, and Insulator #1 in front of Batteries #7 and #8.
4. Disconnect sense leads from J7 and J8 on left side of battery test terminal strip.
5. Disconnect Cable #3 from Battery #7 (+) and slide boot over cable terminal.

6. Disconnect Cable #5 from Battery #8 (–) and slide boot over cable terminal.
7. Disconnect one of the Cable #1s from Battery #7 (–) and Battery #8 (+), and slide a boot over each cable terminal.
8. Remove Batteries #7 and #8 and the separator pad between them.
9. Remove the remaining Insulator #1.

ILLUSTRATION 6-3
BATTERY CONNECTIONS

6-4 Disconnect and Remove Battery Test Terminal Strip

1. Disconnect sense leads from J5 and J6 on left side of battery test terminal strip, and J1, J2, J3 and J4 from right side.
2. Remove battery test terminal strip from battery compartment.

6-5 Continue Battery Disconnection/Removal

1. Disconnect Cable #8 from Battery #5 (+) and slide boot over cable terminal.
2. Disconnect Cable #3 from Battery #6 (–) and slide boot over cable terminal.
3. Disconnect remaining Cable #1 from Battery #6(+) and Battery #5 (–) and slide boot over each cable terminal.
4. Remove Batteries #5 and #6 and the separator pad between them.
5. Remove Insulator #5, #6, and Insulator #2 behind Batteries #5 and #6.
6. Disconnect Cable #4 from Terminal #4 of circuit breaker in AMX1 A3 and slide boot over cable terminal. Disconnect other end of Cable #4 from Battery #1 (+) and slide boot over cable terminal.
7. Disconnect Cable #8 from Battery #4 (–) and slide boot over cable terminal.
8. Disconnect one of the Cable #7s from Battery #1 (–) and Battery #2 (+) and slide a boot over each cable terminal.
9. Disconnect the other Cable #7 from Battery #3 (–) and Battery #4 (+) and slide a boot over each cable terminal.
10. Disconnect Cable #6 from Battery #2 (–) and Battery #3 (+) and slide a boot over each cable terminal.
11. Remove Insulator #4 and the two remaining Insulator #2s.
12. Remove the four remaining batteries and battery separator pads.

6-6 Install Batteries



PERSONAL INJURY ! REMOVE ALL JEWELRY AND RINGS WHEN HANDLING THE BATTERIES OR REACHING IN THE BATTERY COMPARTMENT. USE EXTREME CAUTION, HIGH AMPERAGE AND VOLTAGE IS PRESENT. USE INSULATED TOOLS FOR WORKING ON BATTERIES.

During installation:

- Make all battery connections as shown in Illustration 6–3 using the hardware shown.
- Reference to the right and left sides of battery compartment are relative to the view shown in Illustration 4.
- Cable numbers (#s) indicated in Illustration 6–3 and in this instruction are marked on the parts.
- Insulator numbers are shown for convenience in description, but are not marked on the parts.
- There are no Insulators #3 and #7.
- Battery numbers may or may not be marked on batteries. Mark numbers on batteries if desired.
- Be sure that battery compartment top and bottom insulation is in place.
- Connect cables only in the order instructed, and only when instructed.
- Apply 45–50 pound-inches (5.1–5.6 N-m) torque to ALL battery connections.
- Place insulation boots over ends of each cable.

- Don't remove paper backing from any insulating separator pad. Orient insulator pads with paper backing facing up.

6-7 Prepare to Install Batteries


1. Verify that AMX key switch is OFF, circuit breaker is OFF, and charger cord is unplugged.
2. Ensure that all relevant Service Notes and Field Modification Instructions have been applied.
3. Ensure that new battery voltage is no more than 0.4 volts higher than average battery voltage when replacing individual batteries. For severely discharged sets, 112 Volts or less, recharge the set for one hour before replacing defective battery(s) with new battery(s).

6-8 Install and Connect Batteries

1. Install terminal plates on each battery (see Illustration 6-4). Discard flat washers supplied with batteries.
2. Position Batteries #2 and #4 (see Illustration 6-3) in the battery compartment, one on top of the other with a battery separator pad between them. Put both batteries in back right side of compartment.
3. Position Batteries #1 and #3 in the compartment opening stacked one on top of the other, with a separator pad between them (see Illustration 6-3). Push both batteries into position alongside Batteries #2 and #4.
4. Install Insulator #4 between Batteries #2 and #4 and Batteries #1 and #3.

Note: Don't confuse Insulator #4 with Insulator #1. Insulator #4 has two non-functional holes. Insulator #1s each have two cup-shaped bumps with no holes.

5. Install an Insulator #2 next to Batteries #1 and #3, and another Insulator #2 next to Batteries #2 and #4.

ILLUSTRATION 6-4
INSTALLATION OF TERMINALS AND TERMINAL PLATES

6. Connect end of Cable #6 with sense lead to Battery #3 (+) and connect other end to Battery #2 (-). Slide a boot over each terminal.
7. Connect a Cable #7 end with sense lead to Battery #4 (+), and connect other end to Battery #3 (-). Slide a boot over each terminal.
8. Connect remaining Cable #7 end with sense lead to Battery #2 (+) and connect other end to Battery #1 (-). Slide a boot over each terminal.
9. Connect end of Cable #8 without sense lead to Battery #4 (-), and slide a boot over this terminal. Orient cable so that the terminal lug points down at 45° toward door opening.
10. Connect end of Cable #4 with sense lead to Battery #1 (+) and slide a boot over this terminal. Route Cable #4 through grommet in hole in top of battery compartment, connect it to terminal #4 on circuit breaker in AMX1 A3, and slide boot over terminal.
11. Install Insulator #5, #6, and the remaining Insulator #2.
12. Position Batteries #5 & #6 in compartment opening stacked one on top of the other with a separator pad between them.
13. Connect one of the Cable #1's end with sense lead to Battery #6 (+), and connect other end to Battery #5 (-). Slide boots over each terminal.
14. Connect end of Cable #3 without the sense lead to Battery #6 (-), and slide a boot over terminal.
15. Push Batteries #5 and #6 to rear of battery compartment.
16. Connect free end of Cable #8 (the cable which comes from battery #4 [-]) to Battery #5 (+), and slide a boot over terminal.

6-9 Install and Connect Battery Test Terminal Strip

1. Install the battery test terminal strip (see Illustration 6-2) into battery compartment as shown in Illustration 6-3.
2. Connect Cable #7 Sense Lead (from Battery #2) to Terminal J2 (RHS [right hand side]) as shown in Illustration 6-2.
3. Connect Cable #7 Sense Lead (from Battery #4 [+]) to Terminal J4 (RHS).
4. Connect Cable #6 Sense Lead (from Battery #3 [+]) to J3 (RHS).
5. Connect Cable #4 Sense Lead (from Battery #1 [+]) to J1 (RHS).
6. Connect Cable #1 Sense Lead (from Battery #6 [+]) (LHS [left hand side]).
7. Connect Cable #8 Sense Lead (from Battery #5 [+]) to J5 (LHS).

6-10 Continue Battery Installation/Connection

1. Install the Insulator #1 between Batteries #5 and #6 and Batteries #7 and #8 with cup-shaped bumps pointing away from you. Route free end of Cable #3 over top of Insulator #1.
2. Position Batteries #7 and #8 in compartment opening stacked one on top of the other with a separator pad between them.
3. For remaining Cable #1, connect end with sense lead to Battery #8 (+) and connect other end to Battery #7 (-). Slide a boot over each terminal.
4. Connect end of Cable #5 without sense lead to Battery #8 (-), and slide a boot over terminal.
5. Push Batteries #7 and #8 back into position while holding Cable #3 above Insulator #1.
6. Connect free end of Cable #3 from Battery #6 [-] vertically to Battery #7 (+) and slide a boot over terminal.

6-11 Connect Battery Terminal Test Strip

1. Connect Cable #1 Sense Lead (from Battery #8 [+]) to J8 (LHS).
2. Connect Cable #3 Sense Lead (from Battery #7 [+]) to J7 (LHS).

6-12 Continue Battery Installation/Connection

1. Install another Insulator #1 between Batteries #7 and #8 and Battery #9, with cup-shaped bumps pointing away from you. Route free end of Cable #5 over the top of this Insulator #1.
2. Install Insulator #8.
3. Position Battery #9 on top of Insulator #8.
4. Connect free end of Cable #5 from Battery #8 [-] to Battery #9 (+), and slide a boot over terminal.
5. Tape the end of Cable #2 that does not have the sense lead to prevent it from shorting after it is connected to Battery #9.
6. Connect other end of Cable #2 to Battery #9 (-), and slide a boot over cable terminal. Route Cable #2 through grommet in hole in top of battery compartment.

Note: Do NOT connect any ground cable to the ground stud yet.

6-13 Connect Battery Terminal Test Strip

1. Connect Cable #5 Sense Lead (from Battery #9 [+]) to J9 (LHS).
2. Connect Cable #2 Sense Lead (from Battery #9 [-]) to J10 (LHS).


6-14 Complete Battery Installation/Connection

1. Measure the battery voltage between circuit breaker lead marked #6 and the ground stud. It should measure between 111 and 118 volts DC, depending upon state of charge of the batteries. If it does not, determine the cause of the problem.
2. If unit has a battery sense harness, connect battery sense harness to the battery terminal test strip connector J11. See Illustration 6-2.

6-15 Install Battery Compartment Cover

1. Place the unit on a level floor. Extend the horizontal arm and rotate it to the front end of the AMX (vertical-column end). See Illustration 6-5.

ILLUSTRATION 6-5

BATTERY COMPARTMENT COVER INSTALLATION

Note: Distribution of weight is necessary to install the battery compartment cover.

2. Place the remaining Insulator #1 in the compartment opening with cup-shaped bumps pointing away from you.



The battery compartment cover is an integral part of the Base Assembly. All six bolts must be in place and torqued to 75 pound-inches (8.4 Newton-meters) or the Base Assembly can be warped or cracked.

3. Install battery compartment cover as shown in Illustration 6-5.
4. Remove tape from free ends of all ground cables and connect to ground stud using flat washer and nut.

SECTION 7 – SERVICING THE BASE ASSEMBLY


7-1 Front Bumper

This avoids damage to the front bumper when jacking up the AMX unit.

7-1-1 Remove Pads

1. Remove the 10-32 x 1.25 inch bolts from each pad and remove the corner pads. See Illustration 7-1.
2. Remove the 10-32 x 1.25 inch bolts and flat washers securing the front bumper and remove the middle piece of the front bumper.

ILLUSTRATION 7-1
REMOVAL OF BUMPER ASSEMBLY




7-1-2 Install Pads

1. Install the middle piece of the front bumper using 10-32 x 1.25 inch bolts and flat washers. See Illustration 7-1.
2. Install each of the two corner pads with 10-32 x 1.25 inch bolts for each pad.

7-1-3 Bumper Switch Removal

1. Remove pads. Refer to Section 7-1-1.
2. Remove actuator bracket to release tension on switch actuator. See Illustration 7-2.
3. Remove two binding head screws holding switch on end of bumper track.
4. Remove switch from bumper assembly.

ILLUSTRATION 7-2
BUMPER SWITCH



7-1-4 Bumper Switch Installation

1. Shut off power.
2. Install switch on end of bumper track and secure with binding head screws. See Illustration 7-2.
3. Install actuator bracket to place tension on switch actuator.
4. Verify that bumper switch is activated by spring tension on switch actuator bracket.
5. Install bumper assembly. Refer to Section 6-4.

7-1-5 Bumper Bearings Removal

Note: This procedure covers the removal of any one or more of the eight bearings comprising the bumper assembly.


1. Remove the bumper pads. Refer to Section 7-1-1.
2. Remove the bumper switches. Refer to Section 7-1-3.
3. Remove spring from between actuator bracket and spring bracket to release tension on switch actuator. See Illustration 7-2.
4. Remove switch actuator bracket by removing binding head screws and hex nuts.
5. Slide track from end of bumper assembly to expose bearings.
6. Remove bearing from mounting slot in bumper by removing hex head screw, spacer and shaft.

7-1-6 Bumper Bearings Installation

Note: This procedure covers the removal of any one or more of the eight bearings comprising the bumper assembly.

1. Install bearing in bumper assembly mounting slot and secure with one hex head screw, spacer and shaft. See Illustration 7-3.

ILLUSTRATION 7-3
BUMPER BEARINGS



2. Slide track into end of bumper assembly.
3. Install switch actuator bracket and secure with two binding head screws and hex nuts.
4. Install spring between actuator bracket and spring bracket.
5. Verify that bumper switch is activated by spring tension on switch actuator bracket.
6. Install the bumper assembly. Refer to Section 6-6.

7-2 Raise Unit for Service


Lifting can be done with:

- Automotive floor jack (minimum 1200 lbs. (550 kg) capacity).
- Lever (2 × 4 inches, palette lever) (5 cm × 10 cm).
- Chain hoist.
- AMX-4 Series lifting jack, 46-302966P1.

It is **STRONGLY RECOMMENDED** that the unit be lifted using the AMX-4 Series lifting jack. This instruction is based on the use of this tool.

Get the 46-302966P1 AMX-4 lifting jack and note the two lips on the front of the jack. See Illustration 7-4.

ILLUSTRATION 7-4
AMX-4 LIFTING JACK 46-302966P1



7-2-1 Raise Rear Wheels

1. Rotate front wheel caster toward rear, otherwise unit will be unstable when rear of unit is raised.
2. Block the front wheel and lift drive wheel end of base assembly approximately 10 inches (25 cm) and place it on blocks.
3. Place blocks under wheel being removed to prevent it from falling when axle supports are removed.


7-2-2 Lower Rear Wheels

1. Put the jack under the rear of the AMX unit USING THE TWO LIPS ON THE FRONT OF THE JACK, and lift the AMX unit enough to remove the wood block(s).
2. Lower the AMX unit and remove the jack.

7-2-3 Raise Front Wheels

1. Remove front bumper pads. Refer to Section 7-1-1.
2. Put the jack under the front of the AMX unit. See Illustration 7-5.


ILLUSTRATION 7-5
AMX-4 UNIT WITH JACK IN PLACE (FRONT)



Note: Ensure that the two lips on the front of the jack retain the front of the base in to prevent the unit from falling off the jack.

3. Lift the AMX unit high enough to put the wood block (one 4x4 inch or two 2x4 inch) under the AMX unit as shown in Illustration 7-6. Be sure that the block(s) is (are) centered under the front battery box wall, and be sure that the block(s) extend(s) out at least 2 inches (5 cm) on each side of AMX unit.

ILLUSTRATION 7-6
AMX-4 UNIT SUPPORTED BY WOOD BLOCK (FRONT)



NOTE 1: IF TWO BLOCKS ARE USED, ORIENT BLOCKS FOR MAXIMUM STABILITY.

NOTE 2: CENTER BLOCK(S) WITH FRONT BATTERY BOX WALL.

NOTE 3: LOCATE WOOD BLOCK(S) SO THAT BLOCK(S) STICKS OUT AT LEAST 2 INCHES (5 cm) BEYOND EACH SIDE OF THE AMX UNIT.

7-2-4 Lower Front Wheels

1. Put the jack under the front of the AMX unit USING THE TWO LIPS ON THE FRONT OF THE JACK, and lift the AMX unit enough to remove the wood block(s).
2. Lower the AMX unit and remove the jack.

7-3 Remove Drive Wheels


The following procedure covers removal of either the right or left drive wheel. The wheel mounting support is welded to the side of the base assembly with the notch in the lower surface of the support. All wheels have caliper assemblies.

1. Place the horizontal arm in park position.
2. Remove trim covers. See Section 5.
3. Remove motor and brake assembly. Refer to Section 7-5.



PERSONAL INJURY! WHEN LIFTING THE BACK OF THE UNIT IT CAN BECOME UNSTABLE AND TIP OVER. USE EXTREME CAUTION WHEN PERFORMING THIS OPERATION.

4. Lift the rear of the unit. See Section 7-2.
5. Loosen the four locknuts and flat washers holding the wheel support angle to the base assembly. See Illustration 7-7.
6. Slide the support angle with its wheel block to the side and off of the wheel axle. The mounting holes in the support angle are slotted for this purpose.
7. Leave the drive belt on the axle pulley.
8. Remove the wheel and caliper through the bottom of the base assembly.
9. Remove the four hex socket head capscrews and lockwashers from the top of the wheel mounting support which secure the wheel and caliper assembly to the base assembly.
10. Carefully remove the wheel and caliper assembly through the bottom of the base assembly.

ILLUSTRATION 7-7
DRIVE WHEEL ASSEMBLY

7-4 Install Drive Wheels

Note: **Note:** This procedure covers the installation of either the right or left drive wheel.



PERSONAL INJURY! THE UNIT CAN BE UNSTABLE WHEN ON BLOCKS AND TIP OVER. USE EXTREME CAUTION WHEN PERFORMING THIS OPERATION.

Note: Apply light coat of grease (Lubriplate Fiske 630-A-A) to each end of axle.

1. The unit is on blocks approximately 10 inches (25 cm) above the floor.
2. Place the caliper assembly on the wheel axle.
3. Install the drive belt on the drive wheel pulley.
4. Lift the drive wheel through the bottom of the base assembly and into position in the notch in the bottom of the wheel mounting support. See Illustration 7-7.
5. Attach the wheel and caliper assembly to the wheel mounting support with four hex socket head capscrews. Do not tighten in place.
6. Slide the support angle with its wheel block onto the wheel axle as far as it will go.
7. Tighten the four locknuts and flat washers holding the wheel support angle to the base assembly. Torque to 50–70 pound-inches (5.6 to 7.8 N-m).


8. Secure the wheel and caliper assembly to the wheel mounting support. Torque to 40–60 pound-inches (4.5 to 6.7 N-m).
9. Install splash guard and secure with two lock nuts.
10. Install the motor and brake assembly. Refer to Section 6–10.
11. Lower the unit to the floor.

7-5 Remove Motor and Brake


Note: This procedure covers the removal of either the right or left motor and brake assembly.

1. Move unit to a flat, level surface.
2. Shut off power.
3. Remove cassette drawer from unit. Refer to Section 5–1.
4. Disconnect motor and brake wire connectors from motor control panel. See Illustration 7–8.

ILLUSTRATION 7-8
MOTOR AND BRAKE WIRING



5. Carefully remove spring between base assembly and spring pin on motor mount plate. See Illustration 7–9.
6. Remove nut and bolt holding motor mount to base.
7. Remove drive belt from motor pulley.
8. Remove motor and brake assembly.
9. Remove capscrews and washers holding motor to motor mount plate.
10. Remove binding head screws attaching brake to motor.



7-6 Install Motor and Brake

Note: This procedure covers the installation of either the right or left motor and brake assembly.

1. Shut off power.
2. Install brake on motor and secure with binding head screws. Be certain brake is oriented on motor (wires in position to be connected to motor control panel). See Illustration 7-8.
3. Install motor on mounting plate and attach with hex socket head capscrews and washers.

Note: Motor mounting plate must be oriented with the plate corner bevel toward the rear. If the Spring Pin was removed from the mounting plate, replace it and check that it does not rub on the wheel. See Illustration 7-9.

4. Install motor and brake assembly in the base assembly.
5. Attach motor mount to base with nut and motor mount bolt and tighten finger tight
6. Place drive belt in position on motor pulley.
7. Carefully install spring between base assembly and motor and spring pin on motor mount plate. This places tension on the drive belt.
8. Tighten motor mount nut, ensuring the belt tension is not offset by tightening. Torque to 50–60 foot-pounds (67.8 – 81.3 N-m).

9. Connect motor and brake wire connectors to motor control panel and replace ty-rap.:
J1 – Left Motor J2 – Right Motor


Note: Use caution when installing brake connector. Connector is keyed but can be forced on in reverse direction causing brake and possibly printed wire board damage.

10. Install cassette drawer in frame. Refer to Section 5-2.

7-7 Remove Caster Wheel

1. Remove front bumper. See Section 7-1.
2. Lift front of unit. See Section 7-2.
3. Remove the four nuts and caster. See Illustration 7-10.

ILLUSTRATION 7-10
CASTER REMOVAL/INSTALLATION TYPICAL (RIGHT CASTER SHOWN)



7-8 Install Caster Wheel

1. Install a new caster.
2. Apply Loctite 242 to stud threads and install the four nuts. See Illustration 7-10.
3. Tighten the nuts to 15 to 25 pound-feet (20.5 to 34 N-m).
4. Lower unit. See Section 7-2.
5. Install front bumper. See Section 7-1-2.

SECTION 8 – SERVICING THE DRIVE HANDLE ASSEMBLY


This section provides procedures covering the removal and installation of field serviceable parts in the handle assembly.

8-1 Handle Sensor Removal

Note: This procedure covers removal of either the right or left sensor from the handle assembly.

1. Open top cover and place in position to easily access the handle assembly. Refer to Section 5-5.
2. Disconnect connector at each end of handle assembly. See Illustration 8-1
3. Cut the two ty-raps holding wires to the angle support bracket at each end of the Handle Assembly.
4. Remove handle assembly by removing eight hex flange head capscrews securing it to the two support brackets. See Illustration 8-1.

ILLUSTRATION 8-1
DRIVE HANDLE SUPPORT BRACKET




5. Disconnect sensor leads (red, white and black – positions 3, 4 & 5 in connector) from connector. See Illustration 8-2.

ILLUSTRATION 8-2
DRIVE HANDLE SENSOR LEADS

	1 FROM DRIVE ENABLE SWITCH
	2 FROM DRIVE ENABLE SWITCH
RED	3 FROM SENSOR
WHT	4 FROM SENSOR
BLK	5 FROM SENSOR

WIRING SIDE MALE

6. Detach end cap from end casting by removing four hex socket button head capscrews. See Illustration 8-3.

ILLUSTRATION 8-3
HANDLE ASSEMBLY

7. Slide end cap along handle away from end casting.
8. Remove end casting from handle assembly by removing two hex socket head capscrews securing it to handle springs. Be careful not to damage sensor and drive enable switch leads.
9. Remove two binding head capscrews securing sensor to standoff posts inside the end casting.
10. Pull sensor leads through back of end casting and remove sensor.

8-2 Handle Sensor Installation

Note: This procedure covers installation of either the right or left sensor in the handle assembly.

1. Position sensor on standoff posts inside end casting and secure with two binding head capscrews. See Illustration 8-3.
2. Feed sensor leads through back end of casting to the connector.
3. Position handle assembly in end casting and secure handle springs to casting with two hex socket head capscrews. Be careful not to damage sensor and drive enable switch leads.


4. Slide end cap along handle and position on end casting.
5. Secure end cap on end casting with four hex socket button head capscrews.
6. Connect sensor leads (red, white and black – positions 3, 4 & 5 in connector) to connector). See Illustration 8-2.
7. Fasten handle assembly to two support brackets with eight hex flange head capscrews. See Illustration 8-1.
8. Connect connectors at each end of handle assembly.
9. Connect ground wire, located at left end of handle, to mainframe.
10. Attach wires to the handle assembly angle support bracket with ty-raps.
11. Close top cover. Refer to Section 5-6.

8-3 Handle Drive Enable Switch Removal

Note: This procedure covers removal of either the right or left enable switch from the handle assembly.

1. Open top cover and place in position to easily access the handle assembly. Refer to Section 5-5.
2. Cut ty-raps holding drive enable switch leads to end casting angle support bracket. See Illustration 8-1.
3. Free the drive enable switch leads from the other wires. (blue, yellow or gray – positions 1 & 2 in connector) See Illustration 8-2.
4. Disconnect drive enable switch leads from connector.
5. Remove end cap from end casting by removing four hex socket button head capscrews. See Illustration 8-3.
6. Slide end cap along handle away from end casting to allow access to drive enable switch.
7. Pull switch leads through back of end casting.
8. See Illustration 8-4 for the drive enable microswitch.
9. Remove two screws securing drive enable microswitch to bracket, and remove this microswitch.

ILLUSTRATION 8-4
DRIVE ENABLE MICROSWITCH



8-4 Handle Drive Enable Switch Installation

Note: This procedure covers installation of either the right or left drive enable switch in the handle assembly.

1. See Illustration 8-3, Detail "B", for the drive enable microswitch.
2. Place the drive enable microswitch on its bracket and secure it in place with two screws.
3. Feed drive enable switch leads through back of end casting to connector.

4. Slide end cap along handle into position on end casting.
5. Secure end cap to end casting with four hex socket button head capscrews.
6. Connect drive enable switch leads (blue, yellow or gray) to connector (positions 1 & 2 in connector). See Illustration 8-2.
7. Secure enable switch leads and other leads to end casting support bracket with a ty-rap. See Illustration 8-1.
8. Close top cover. Refer to Section 5-6.


SECTION 9 – SERVICING THE HIGH VOLTAGE TRANSFORMER

This section provides procedures covering the removal and installation of the high voltage transformer.

9-1 High Voltage Transformer Removal

1. Shut power off by turning the circuit breaker OFF.
2. Remove front cover by carefully pulling cover from four ball stud retainer sockets.
3. Disconnect anode and cathode cables from high voltage transformer.
4. Remove two cap screws holding rear cover on base assembly and remove cover.
5. Remove cassette drawer. Refer to Section 5-1.
6. Remove cassette drawer magnetic strip lock by removing three binding head screws and flat washers.
7. Remove high voltage transformer terminal cover by removing two binding head screws and washers. See Illustration 9-1.

ILLUSTRATION 9-1
HIGH VOLTAGE TRANSFORMER



8. Tag and label all leads on terminal
9. Disconnect all leads from terminal.
10. Tag and label all connections to Drive Electronics Module Assembly.
11. Disconnect all connections to Drive Electronics Module Assembly.
12. Remove Drive Electronics Assembly by removing four hex washer head screws securing it to mainframe.



PERSONAL INJURY! HIGH VOLTAGE TRANSFORMER WEIGHS APPROXIMATELY 100 LBS (45KG). USE CAUTION WHEN MOVING IT.

Note: Before sliding high voltage transformer to rear of unit, move cable harness out of way of transformer cap.

13. Use two people to remove high voltage transformer.
14. While one person pushes transformer from front of unit, second person guides transformer out back of unit until it rests on rear edge of base assembly.
15. Use two people to lift transformer out of unit.

9-2 High Voltage Transformer Installation

1. Shut power off by turning circuit breaker OFF.
2. Use two people to install transformer.



PERSONAL INJURY! HIGH VOLTAGE TRANSFORMER WEIGHS APPROXIMATELY 100 LBS (45KG). USE CAUTION WHEN MOVING IT.

3. Lift transformer and place front side on edge of opening into mainframe. Rest the terminal side of transformer on rear edge of base assembly.

Note: Before sliding high voltage transformer into unit:

- Check oil level
 - Check vent is open
 - Set spark gap to 0.01 inches (.25 millimeters)
 - Move cable harness out of way of transformer cap.
4. While one person pushes transformer from back of unit, second person guides transformer into unit until it contacts transformer restraint plate.
 5. Install Drive Electronics Module Assembly and secure it to mainframe with four hex washer head screws.
 6. Connect all connections to Drive Electronics Module Assembly in accordance with tagged and labeled instructions on connectors.
 7. Connect all leads to transformer terminal in accordance with tagged and labeled instructions on leads.
 8. Remove all tags.
 9. Install high voltage transformer terminal cover and secure in place with two binding head screws and washers.
 10. Install cassette drawer magnetic strip lock and secure with three binding head screws and flat washers.
 11. Install cassette drawer. Refer to Section 5-2.
 12. Install rear cover and secure with two flat head cap screws.
 13. Connect anode and cathode cables to high voltage transformer.
 14. Install front cover with two flat head capscrews.

SECTION 10 – SERVICING THE POWER CORD AND REEL ASSEMBLY AND CABLE ROUTING


This section covers the removal and installation of the main power cord reel, and the horizontal arm/vertical column cable routing.

10-1 Cord/Cord Reel Removal

10-1-1 Cord Reel Removal (Domestic Models)

1. Shut power off by turning circuit breaker OFF.
2. Open top cover. Refer to Section 5-5.
3. Remove four hex nuts attaching low voltage power supply to mainframe and slide it backward out of mainframe housing. See Illustration 10-1.
4. Disconnect cord reel leads from F1, F2 and ground stud.
5. Remove two internal hex cap screws retaining grommet in cut-out in mainframe and remove cable and grommet.
6. Remove two hex washer head screws and lock nuts holding cord reel to side of low voltage power supply housing.
7. Lift off cord reel pulling leads through grommet in low voltage power supply housing.

ILLUSTRATION 10-1
CORD REEL REMOVAL/INSTALLATION




10-1-2 Cord Reel Removal (IEC Models)

1. Shut power off by turning circuit breaker OFF.
2. Open top cover. Refer to Section 5-5.
3. Disconnect cord reel leads at main terminal block. See Illustration 10-2.
4. Remove two internal hex cap screws retaining grommet in cut-out in mainframe and remove cable and grommet.
5. Remove two hex washer head screws and lock nuts holding cord reel to side of low voltage power supply housing.
6. Loosen screw in strain relief bracket located behind cord reel.

7. Lift off cord reel pulling leads through strain relief bracket.

ILLUSTRATION 10-2
CORD REEL REMOVAL/INSTALLATION



10-2 Cord/Cord Reel Installation

10-2-1 Cord Reel Installation (Domestic Models)

1. Shut power off by turning circuit breaker OFF.
2. Feed cord reel leads through grommet in low voltage power supply housing. See Illustration 10-1.
3. Install cord reel on side of low voltage power supply housing and secure with two hex washer head screws and lock nuts.
4. Install cable grommet in cut-out in mainframe with two cap screws.
5. Connect three leads from cord reel to F1 (Black), F2 (White) and ground stud (Green).
6. Slide low voltage power supply into mainframe housing and fasten with four hex nuts.
7. Close top cover. Refer to Section 5-6.

10-2-2 Cord Reel Installation (IEC Models)

1. Shut power off by turning circuit breaker OFF.
2. Feed cord reel leads through strain relief bracket and tighten bracket retaining screw. See Illustration 10-2.
3. Install cord reel on side of low voltage power supply housing and secure with two hex washer head screws and lock nuts.
4. Install cable grommet in cut-out in mainframe with two cap screws.
5. Connect three leads from cord reel to main terminal block. Use color code as identified by tag adjacent to terminal block.
6. Close top cover. Refer to Section 5-6.

10-3 Horizontal Arm and Vertical Column Cable Routing

There are seven cables mounted on the horizontal arm and column. These must be carefully routed and mounted in order for the equipment to be operated with maximum flexibility and to prevent binding or jamming of cables which would restrict positioning the X-ray tube.

Cables:


- Cathode – X-ray tube to transformer located in base of unit.
- Anode – X-ray tube to transformer located in base of unit.
- Collimator – Collimator to terminal box located on vertical column assembly
- X-ray Control – X-ray tube to terminal box
- Horizontal Arm Brake – Horizontal arm to terminal box
- Control – Terminal box to base unit
- Upper Column Brake – Upper brake to terminal box

When connecting the anode or cathode cables to the tube or transformer, see Direction 2196272-100, *High Voltage Cable Installation and Troubleshooting Procedures*.

10-3-1 Cathode Cable


1. The cathode cable starts at the tube cathode connector, runs through the yoke, crosses underneath the opposite side of the yoke to clamp point "K." See Illustration 10-3.

ILLUSTRATION 10-3

CATHODE CABLE ROUTING FROM TUBE TO FIRST CLAMP


2. From clamp "K" the cable continues to strap point "M" and then to clamp point "H." See Illustration 10-4.

ILLUSTRATION 10-4
CATHODE CABLE ROUTING FROM TUBE



3. The cathode cable continues from clamp point "H" and strap point "G," to clamp point "E" to the front panel in a figure-eight configuration, and then terminates at the transformer in the base of the unit.
- > The distances are referenced in Illustration 10-5.

ILLUSTRATION 10-5
CATHODE CABLE




Note: Strap point "J" is installed with collimator cable.

10-3-2 Anode Cable


1. The anode cable starts at the tube anode connector, runs through the yoke crosses underneath to the opposite side of the yoke to clamp point "D." See Illustration 10-6.

ILLUSTRATION 10-6
ANODE CABLE ROUTING FROM TUBE




2. From clamp point "L" the cable continues to strap point "N" and clamp point "D." See Illustration 10-7.

ILLUSTRATION 10-7
ANODE CABLE ROUTING FROM FIRST CLAMP TO TRANSFORMER



- > The distance from the tube to each clamp and strap point is shown in Illustration 10-8.

ILLUSTRATION 10-8
ANODE CABLE




3. The anode cable continues from clamp point "D" to strap point "C," to point "A" and then to the front panel in a figure-eight configuration, then terminates at the transformer in the base of the unit.

10-3-3 Collimator Cable


1. The collimator cable starts at the collimator, passes over the anode cable and goes to strap point "J" where it is joined to the cathode cable. There is a large loop of cable before reaching strap point "J" as shown in Illustration 10-9.

ILLUSTRATION 10-9
COLLIMATOR CABLE ROUTE TO FIRST STRAP




2. Continue routing and ty-rapping the collimator cable along the cathode cable to strap point "C." See Illustration 10-10.

ILLUSTRATION 10-10
CATHODE AND COLLIMATOR CABLES




3. There is a large loop of collimator cable immediately after strap point "C." See Illustration 10-11.

ILLUSTRATION 10-11
COLLIMATOR CABLE ROUTE ON CATHODE CABLE



4. The collimator cable goes from strap point "C" directly into the vertical column assembly terminal box where it is connected to terminal strip AMX3 A1 TS1. See Illustration 10-12 and Table 10-1.

ILLUSTRATION 10-12
COLUMN TERMINAL STRIP AMX3 A1 TS1TABLE 10-1
AMX3 A1 TS1 CONNECTIONS


Terminal No.	Signal Name	FROM		TO		
		AMX1 A2 TS3	Color	Color	Terminal	
1	ROTOR PHASE2	3	BLK	BLK	AMX3 A3 TS1-5	X-RAY TUBE
2	ROTOR COMN	4	WHT	WHT	AMX3 A3 TS1-3	X-RAY TUBE
3	ROTOR PHASE 1	5	GRN	GRN	AMX3 A3 TS1-4	X-RAY TUBE
4	X-RAY PRES SW	6	WHT/ORN	BLU	AMX3 A3 TS1-1	X-RAY TUBE
5	X-RAY PRES SW RTN	7	GRY	RED	AMX3 A3 TS1-2	X-RAY TUBE
6	FIELD LIGHT SW	8	WHT/BRN	WHT	AMX3 A4	COLLIMATOR
7	FIELD LIGHT SW RTN	9	WHT/RED	YEL	AMX3 A4	COLLIMATOR
8	+12VDC LOCKS	10	WHT/BLK	RED	AMX3 BRK2	COLUMN LOCK
9	– JUMPER FROM 8 –			RED	AMX3 A2 J1	TUBE ARM LOCK
10	–BATT	11	VIO	BLK	AMX3 BRK2	COLUMN LOCK RTN
11	– JUMPER FROM 10 –			BLK	AMX3 A2 J1	TUBE ARM LOCK RTN
12	PWR GND	12	YEL	ORN	AMX3 A4	COLLIMATOR
13	+24VDC	13	BLU	GRN	AMX3 A4	COLLIMATOR
14	FIELD LIGHT	14	BRN	RED/BLK	AMX3 A4	COLLIMATOR
15	FIELD LIGHT RTN	15	ORN	BRN/BLU	AMX3 A4	COLLIMATOR
16	NOT USED					

Terminal No.	Signal Name	AMX1 A2 TS3	Color	Color	Terminal	
GND STUD		BASE GND STUD	GRN/YEL	GRN/YEL	AMX3 A2 J1	TUBE ARM GND
				GRN/YEL	AMX3 A4	COLLIMATOR GND

10-3-4 Stator Cable


1. The stator cable starts at the tube, passes under the cathode cable and joins the anode cable just before clamp "L." There is a large loop of cable before reaching clamp point "L." The cables are ty-rapped on each side of clamp point "L" as shown in Illustration 10-13.

ILLUSTRATION 10-13
STATOR CABLE ROUTE TO FIRST TY-RAP




2. The stator cable follows the anode cable from clamp point "L" to strap point "M" with two ty-raps in between. See Illustration 10-14.

ILLUSTRATION 10-14
STATOR CABLE ROUTE ON ANODE CABLE



3. Continue routing and ty-rapping the collimator cable along the anode cable to strap point "G." See Illustration 10-15.

ILLUSTRATION 10-15
ANODE AND STATOR CABLES



4. The horizontal arm brake cable is also picked up and retained until strap point "G." See Illustration 10-14.


5. There is a large loop after strap point “G” before the cable enters the bottom of the vertical column assembly terminal box and is connected to the terminal strip AMX# A1 TS1. See Illustration 10-12 and Table 10-1.

10-3-5 Horizontal Arm Brake Cable

1. The horizontal arm brake cable exits the left side of the horizontal arm outer tube and is ty-rapped to the bottom of the anode cable just below clamp point “H.” See Illustration 10-16.

ILLUSTRATION 10-16


HORIZONTAL ARM BRAKE ROUTING



2. The cable continues from near clamp point “H” and is ty-rapped to the anode and stator cables at the locations shown in Illustration 10-16 until reaching strap point “G”.
3. There is a large loop after strap point “G” before the cable enters the bottom of the vertical column assembly terminal box and is connected to the terminal strip AMX3 A1 TS1. See Illustration 10-12 and Table 10-1.


10-3-6 Control Cable

The control cable exits the bottom of the terminal box on the rear of the vertical column assembly, passes beneath the anode cable and connects directly into the unit just below the anode and cathode cables. The cable is loosely draped to accommodate arm and column movements. See Illustration 10-17.

ILLUSTRATION 10-17
CONTROL CABLE ROUTING

10-3-7 Vertical Column Brake Cable

1. The vertical column brake cable exits the side of the vertical column assembly and is routed to a channel on the rear of the column. See Illustration 10-18.

ILLUSTRATION 10-18
UPPER BRAKE CABLE ROUTE ON COLUMN

2. The brake cable enters the channel at the top of the vertical column assembly and continues down into the top of the terminal box and is connected to the terminal strip AMX3 A1 TS1. See Illustration 10-12 and Table 10-1.



GE Medical Systems

*GE Medical Systems: Telex 3797371
P.O. Box 414, Milwaukee, Wisconsin 53201 U.S.A.
(Asia, Pacific, Latin America, North America)*

*GE Medical Systems — Europe: Telex 698626
283, rue de la Minière, B.P. 34, 78533 Buc Cedex
France*